

DATASHEET

This Tracheotomy Trainer is usable for various methods of tracheotomy.

It allows training of both conventional and percutaneous tracheotomy and puncture/incision into the cricothyroid ligament.

For conventional tracheotomy, different types of incisions can be performed including longitudinal, transverse, crucial and (inverted) U-shaped.

The trainer allows user to determine the proper incision site while identifying arterial locations, and to observe intratracheal conditions from the head.

Reusable trachea and skin

The trachea can be reused by turning it upside down or inside out, and the skin by sliding it to either side or turning it upside down

Tracheotomy Trainer



KEY BENEFITS

- Training of both conventional and percutaneous tracheotomy and puncture/incision into the cricothyroid ligament
- For conventional tracheotomy: longitudinal, transverse, crucial, (inverted) U-shaped incisions possible
- Determination of proper incision site and observation of intratracheal conditions
- Reusable trachea and skin

ORDERING SPECIFICATIONS

Item no.

Tracheotomy Trainer

SB45067U incl. torso, neck with arteries, thyroid cartilage with cricoid cartilage, 4 x trachea, 4 x skin, lubricant

SB50107U Skin Tracheotomy Trainer Set 3

SB50108U Trachea for Tracheotomy Trainer Set 3

Ambu

■ Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com