

DATASHEET

The Complete GERi is a geriatric patient care simulator with natural flexibility and superior range of motion.

It has an overall female appearance but is also convertible to a male with removal of the wig and attachment of male genitals.

The simulator enables the practice of bandaging and wound dressing, bed baths and hair care washing, clothing changes, or ear canal irrigation with otic drops administration and hearing aid placement as well as several medical and nursing procedures.

These are the training of oral/nasal hygiene like lavage, gavage and suctioning, the placement of a nasogastric tube, and intramuscular injections on arm, thigh and buttock.

For visual inspection the manikin features normal and cancerous mole, stage 1 sacral ulcer, dilated and constricted pupils and reddened skin folds. The GERi manikin allows also training of ostomy and tracheostomy care, as well as gastrostomy procedures, lavage and gavage or enema administration (female only).

The urinary catheterization can be practised on female and uncircumcised male.

The male genitalia features also prostate exam.

Complete GERi™



KEY BENEFITS

- Training manikin that is helpful for preparing personnel for child rescue where highly charged emotional factors come into play
- Patient care simulator with elderly appearance, skin wrinkles and folds
- All major joints move natural for practicing patient positioning with realistic flexibility and range of motion
- Simulating over 35 nursing and medical procedures
- Internal fluid reservoirs
- Tracheostomy care
- Eye / Ear-canal irrigation treatment
- Intramuscular injection sites
- NG tube placement
- Ostomy care, lavage and suctioning
- Gastrostomy procedures
- Urinary catheterization

ORDERING SPECIFICATIONS

Item no.

Complete GERi

LF04001

Fullbody manikin incl. male/female genitalia, wig, dentures, hearing aid, syringes, lubricant, fluid drainage basin, stomach reservoir bag, bladder reservoir bag, bladder reservoir pressure sleeve, catheter, enema reservoir bag

LF04080

IV Arm Replaceable skins and latex veins that actually roll when palpated

LF04079

Blood Pressure Arm for GERi™ Reproduces the five Korotkoff sounds and allows the instructor to vary systolic and diastolic levels, pulse rate, volume, and auscultatory gap.

Ambu

Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com