




PUB-000227
Symbol Explanation
ISO 15223-1




Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
REF	5.1.6	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the manufacturer's <i>catalogue number</i> so that the <i>medical device</i> can be identified	English	Catalogue number
				Bulgarian	Каталожен номер
				Croatian	Kataloški broj
				Chinese	目录号
				Czech	Katalogové číslo
				Danish	Katalognummer
				Dutch	Catalogusnummer
				Estonian	Katalooginumber
				Finnish	Luettelonumero
				French	Référence catalogue
				German	Katalognummer
				Greek	Αριθμός καταλόγου
				Hungarian	Katalógusszám
				Italian	Numero di catalogo
				Japanese	カタログ番号
				Latvian	Kataloga numurs
				Lithuanian	Katalogo numeris
				Norwegian	Katalognummer
				Polish	Numer katalogowy
				Portuguese	Número de catálogo
Romanian	Număr de catalog				
Slovak	Katalógové číslo				
Slovenian	Kataloška številka				
Spanish	Número de catálogo				
Swedish	Katalognummer				
Turkish	Katalog numarası				



Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
LOT	5.1.5	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the manufacturer's batch code so that the batch or lot can be identified	English	Batch code
				Bulgarian	Партиден код
				Croatian	Kod serije
				Chinese	批代码
				Czech	Kód šarže
				Danish	Batchkode
				Dutch	Lotnummer
				Estonian	Partitähis
				Finnish	Eräkoodi
				French	Code du lot
				German	Chargennummer
				Greek	Κωδικός παρτίδας
				Hungarian	Gyártási tétel kódja
				Italian	Codice lotto
				Japanese	バッチコード
				Latvian	Partijas kods
				Lithuanian	Partijos kodas
				Norwegian	Partikode
				Polish	Kod partii
				Portuguese	Código do lote
Romanian	Codul lotului				
Slovak	Kód šarže				
Slovenian	Koda proizvodne serije				
Spanish	Código de lote				
Swedish	Batchnummer				
Turkish	Parti kodu				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.1.1	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the <i>medical device manufacturer</i>	English	Manufacturer
				Bulgarian	Производител
				Croatian	Proizvođač
				Chinese	制造商
				Czech	Výrobce
				Danish	Producent
				Dutch	Fabrikant
				Estonian	Tootja
				Finnish	Valmistaja
				French	Fabricant
				German	Hersteller
				Greek	Κατασκευαστής
				Hungarian	Gyártó
				Italian	Produttore
				Japanese	製造者
				Latvian	Ražotājs
				Lithuanian	Gamintojas
				Norwegian	Produsent
				Polish	Producent
				Portuguese	Fabricante
Romanian	Fabricant				
Slovak	Výrobca				
Slovenian	Proizvajalec				
Spanish	Fabricante				
Swedish	Tillverkare				
Turkish	Üretici				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.1.3	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the date when the <i>medical device</i> was manufactured	English	Date of manufacture
				Bulgarian	Дата на производство
				Croatian	Datum proizvodnje
				Chinese	制造日期
				Czech	Datum výroby
				Danish	Fremstillingsdato
				Dutch	Productiedatum
				Estonian	Tootmise kuupäev.
				Finnish	Valmistuspäivä
				French	Date de fabrication
				German	Herstellungsdatum
				Greek	Ημερομηνία κατασκευής
				Hungarian	A gyártás ideje
				Italian	Data di produzione
				Japanese	製造日付
				Latvian	Ražošanas datums
				Lithuanian	Gamybos data
				Norwegian	Produksjonsdato
				Polish	Data produkcji
				Portuguese	Data de fabrico
Romanian	Data fabricației				
Slovak	Dátum výroby				
Slovenian	Datum proizvodnje				
Spanish	Fecha de fabricación				
Swedish	Tillverkningsdatum				
Turkish	Üretim tarihi				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.1.4	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the date after which the <i>medical device</i> is not to be used	English	Use-by date
				Bulgarian	Срок на годност
				Croatian	Rok trajanja
				Chinese	有效期
				Czech	Datum použitelnosti
				Danish	Sidste anvendelsesdato
				Dutch	Uiterste gebruiksdatum
				Estonian	Aegumiskuupäev
				Finnish	Viimeinen käyttöpäivä
				French	Date de péremption
				German	Haltbarkeitsdatum
				Greek	Ημερομηνία λήξης
				Hungarian	Felhasználható
				Italian	Utilizzare entro (data)
				Japanese	使用期限
				Latvian	Izlietot līdz
				Lithuanian	Naudojimo terminas
				Norwegian	Må brukes innen dato
				Polish	Data ważności
				Portuguese	Data de validade
Romanian	Data de expirare				
Slovak	Dátum spotreby				
Slovenian	Rok trajanja				
Spanish	Fecha de caducidad				
Swedish	Bäst-före-datum				
Turkish	Son kullanma tarihi				



Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 <p>Examples</p>  	5.1.11	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	To identify the country of manufacture of products. Ambu also uses this symbol to state the date of manufacture	English	Country of manufacture. Date of manufacture
				Bulgarian	Държава на производство, Дата на производство
				Croatian	Zemlja proizvodnje, Datum proizvodnje
				Chinese	原产国, 生产日期
				Czech	Země výrobcce, Datum výroby
				Danish	Producentland, Fremstillingsdato
				Dutch	Land van fabrikant, Productiedatum
				Estonian	Tootja riik, Tootmise kuupäev
				Finnish	Valmistusmaa, Valmistuspäivä
				French	Pays de fabrication, Date de fabrication
				German	Herstellungsland, Herstellungsdatum
				Greek	Χώρα κατασκευής, Ημερομηνία κατασκευής
				Hungarian	A gyártó országa, A gyártás ideje
				Italian	Paese di produzione, Data di produzione
				Japanese	製造国, 製造日付
				Latvian	Ražotājvalsts, Ražošanas datums
				Lithuanian	Gamintojo šalis, Gamybos data
				Norwegian	Produksjonsland, Produksjonsdato
				Polish	Kraj produkcji, Data produkcji
				Portuguese	País de fabrico, Data de fabrico
Romanian	Țara de fabricație, Data fabricației				
Slovak	Krajina výroby, Dátum výroby				
Slovenian	Država proizvodnje, Datum proizvodnje				
Spanish	País de fabricación, Fecha de fabricación				
Swedish	Tillverkningsland, Tillverkningsdatum				
Turkish	Üretildiği ülke, Üretim tarihi				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 <p>Example:</p> 	5.3.7	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the temperature limits to which the <i>medical device</i> can be safely exposed	English	Temperature limit
				Bulgarian	Температурно ограничение
				Croatian	Temperaturna ograničenja
				Chinese	温度限制
				Czech	Teplotní limit
				Danish	Temperaturgrænse
				Dutch	Temperatuurlimiet
				Estonian	Temperatuuri piirang
				Finnish	Lämpötilan raja-arvot
				French	Limite de température
				German	Temperaturbegrenzung
				Greek	Όριο θερμοκρασίας
				Hungarian	Hőmérsékleti korlátozás
				Italian	Limite temperatura
				Japanese	温度制限
				Latvian	Temperatūras ierobežojums
				Lithuanian	Temperatūros ribos
				Norwegian	Temperaturgrense
				Polish	Zakres temperatury
				Portuguese	Limite de temperatura
Romanian	Limite de temperatură				
Slovak	Teplotný limit				
Slovenian	Omejitev temperature				
Spanish	Límite de temperatura				
Swedish	Temperaturgräns				
Turkish	Sıcaklık sınırı				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.4.2	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a <i>medical device</i> that is for one <i>single use only</i>	English	Do not re-use
				Bulgarian	Не използвайте повторно
				Croatian	Nemojte ponovno upotrebljavati
				Chinese	请勿重复使用
				Czech	Nepoužívejte opakovaně
				Danish	Må ikke genbruges
				Dutch	Niet opnieuw gebruiken
				Estonian	Mitte korduskasutada / Ainult ühekordseks kasutamiseks
				Finnish	Kertakäyttöinen
				French	Ne pas réutiliser
				German	Nur für den einmaligen Gebrauch
				Greek	Μην επαναχρησιμοποιείτε
				Hungarian	Ne használja újra!
				Italian	Non riutilizzare
				Japanese	再使用しないでください
				Latvian	Nelietot atkārtoti
				Lithuanian	Nenaudoti pakartotinai
				Norwegian	Skal ikke gjenbrukes
				Polish	Nie używać ponownie
				Portuguese	Não reutilizar
Romanian	A nu se reutiliza				
Slovak	Nepoužívajte opakovane				
Slovenian	Za enkratno uporabo				
Spanish	No reutilizar				
Swedish	Får inte återanvändas				
Turkish	Tekrar kullanmayın				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.4.12	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a <i>medical device</i> that may be used multiple times (multiple procedures) on a single patient	English	Single patient multiple use
				Bulgarian	За многократна употреба от един пациент
				Croatian	Za višekratnu upotrebu na jednom pacijentu
				Chinese	可供单个患者多次使用
				Czech	K opakovanému použití u jednoho pacienta
				Danish	Flergangsbrug til én enkelt patient
				Dutch	Meervoudig gebruik bij één patiënt
				Estonian	Korduskasutuseks ühel patsiendil
				Finnish	Moninkertaiseen käyttöön yhdellä potilaalla
				French	Peut être utilisé plusieurs fois surchez le même patient
				German	Einzelner Patient – mehrfach verwendbar
				Greek	Για πολλαπλές χρήσεις σε έναν ασθενή
				Hungarian	Egyetlen pácienshez, több alkalommal használható eszköz
				Italian	Monopaziente, multiuso
				Japanese	1人の患者が複数回使用できる
				Latvian	Atkārtotai lietošanai vienam pacientam
				Lithuanian	Daugkartinio naudojimo vienam pacientui
				Norwegian	Til bruk flere ganger for én pasient
				Polish	Produkt wielokrotnego użytku, dla jednego pacjenta
				Portuguese	Utilização múltipla num único paciente
Romanian	Utilizare multiplă la un singur pacient				
Slovak	Na opakované použitie u jedného pacienta				
Slovenian	Za večkratno uporabo pri enem bolniku				
Spanish	Varios usos en un solo paciente				
Swedish	För flergångsbruk på en patient				
Turkish	Tek hastada tekrar kullanılabilir				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
STERILE R	5.2.4	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a <i>medical device</i> that has been sterilized using irradiation	English	Sterilized using irradiation
				Bulgarian	Стерилизирано чрез облъчване
				Croatian	Sterilizirano zračenjem
				Chinese	采用辐照灭菌
				Czech	Sterilizováno zářením
				Danish	Strålesteriliseret
				Dutch	Gesteriliseerd door bestraling
				Estonian	Steriliseeritud kiirgusega
				Finnish	Steriloitu säteilyttämällä
				French	Stérilisé par irradiation
				German	Sterilisiert durch Bestrahlung
				Greek	Αποστειρωμένο με χρήση ακτινοβολίας
				Hungarian	Besugárzással sterilizálva
				Italian	Sterilizzato con irradiazione
				Japanese	放射線滅菌
				Latvian	Sterilizēts, izmantojot apstarošanu
				Lithuanian	Sterilizuota spinduliuote
				Norwegian	Sterilisert med stråling
				Polish	Steryliзовany radiacyjnie
				Portuguese	Esterilizado por irradiação
Romanian	Sterilizat prin iradiere				
Slovak	Sterilizované ožarovaním				
Slovenian	Sterilizirano z obsevanjem				
Spanish	Esterilización con irradiación				
Swedish	Steriliserad med strålning				
Turkish	Radyasyonla sterilize edilmiştir				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 <p>Example: in combination with 5.2.4</p> 	5.2.11	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a single <i>sterile</i> barrier system	English	Single sterile barrier system
				Bulgarian	Единична стерилна преградна система
				Croatian	Sustav jednostruke sterilne barijere
				Chinese	单层无菌屏障系统
				Czech	System jedné sterilní bariéry
				Danish	Enkelt sterilt barrieresystem
				Dutch	Enkelvoudig steriel barrièresysteem
				Estonian	Ühekordne steriilne barjäärisüsteem
				Finnish	Yksinkertainen steriili estejärjestelmä
				French	Système de barrière stérile unique
				German	System mit nur einer Sterilbarriere
				Greek	Σύστημα μονού στείρου φραγμού
				Hungarian	Egyszeres sterilgát-rendszer
				Italian	Sistema a barriera sterile singola
				Japanese	シングル滅菌バリアシステム
				Latvian	Vienas sterilas barjeras sistēma
				Lithuanian	Viengubo sterilaus barjero sistema
				Norwegian	Enkelt, steril barrieresystem
				Polish	System pojedynczej osłony sterylnej
				Portuguese	Sistema de barreira estéril único
				Romanian	Sistem cu barieră sterilă unică
Slovak	System jednej sterilnej bariéry				
Slovenian	Enojni sterilni pregradni sistem				
Spanish	Sistema de barrera estéril única				
Swedish	Enkelt sterilbarriärsystem				
Turkish	Tekli steril bariyer sistemi				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
<div style="border: 1px solid black; display: inline-block; padding: 2px;">STERILE</div> <div style="border: 1px solid black; display: inline-block; padding: 2px; margin-left: 5px;">EO</div>	5.2.3	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a <i>medical device</i> that has been sterilized using ethylene oxide	English	Sterilized using ethylene oxide
				Bulgarian	Стерилизирано с етиленов оксид
				Croatian	Sterilizirano etilen-oksidadom
				Chinese	使用环氧乙烷杀菌
				Czech	Sterilizováno etylenoxidem
				Danish	Steriliseret med etylenoxid
				Dutch	Gesteriliseerd met behulp van ethyleenoxide
				Estonian	Steriliseeritud etüleenoksiidiga
				Finnish	Steriloitu etyleenioksidilla
				French	Stérilisé à l'oxyde d'éthylène
				German	Mit Ethylenoxid sterilisiert
				Greek	Αποστειρωμένο με τη χρήση αιθυλενοξειδίου
				Hungarian	Etilén-oxidál sterilizálva
				Italian	Sterilizzato con ossido di etilene
				Japanese	エチレンオキサイドによる滅菌済み
				Latvian	Sterilizēts ar etilēna oksīdu
				Lithuanian	Sterilizuota naudojant etileno oksidą
				Norwegian	Sterilisert med etylenoksid
				Polish	Steryliwowana przy użyciu tlenku etylenu
				Portuguese	Esterilizado com óxido de etileno
Romanian	Sterilizat cu etilenoxid				
Slovak	Sterilizované etylénoxidom				
Slovenian	Sterilizirano z etilen oksidom				
Spanish	Esterilización con óxido de etileno				
Swedish	Steriliserad med etylenoxid				
Turkish	Etilen oksitle sterilize edilmiştir				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.2.11 in combination with 5.2.3	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a single <i>sterile</i> barrier system	English	Single sterile barrier system
				Bulgarian	Единична стерилна преградна система
				Croatian	Sustav jednostruke sterilne barijere
				Chinese	单层无菌屏障系统
				Czech	Systém jedné sterilní bariéry
				Danish	Enkelt sterilt barrieresystem
				Dutch	Enkelvoudig steriel barrièresysteem
				Estonian	Ühekordne steriilne barjäärisüsteem
				Finnish	Yksinkertainen steriili estejärjestelmä
				French	Système de barrière stérile unique
				German	System mit nur einer Sterilbarriere
				Greek	Σύστημα μονού στείρου φραγμού
				Hungarian	Egyszeres sterilgát-rendszer
				Italian	Sistema a barriera sterile singola
				Japanese	エチレンオキサイドによる滅菌済み
				Latvian	Vienas sterilas barjeras sistēma
				Lithuanian	Viengubo sterilaus barjero sistema
				Norwegian	Enkelt, sterilt barrieresystem
				Polish	System pojedynczej osłony sterylnej
				Portuguese	Sistema de barreira estéril único
Romanian	Sistem cu barieră sterilă unică				
Slovak	Systém jednej sterilnej bariéry				
Slovenian	Enojni sterilni pregradni sistem				
Spanish	Sistema de barrera estéril única				
Swedish	Enkelt sterilbarriärsystem				
Turkish	Tekli steril bariyer sistemi				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.2.7	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a medical device that has not been subjected to a sterilization process	English	Non-sterile
				Bulgarian	Нестерилно
				Croatian	Nesterilno
				Chinese	未消毒
				Czech	Nesterilní
				Danish	Usteril
				Dutch	Niet-steriel
				Estonian	Mittesteriilne
				Finnish	Steriloimaton
				French	Non stérile
				German	Unsteril
				Greek	Μη αποστειρωμένο
				Hungarian	Nem steril
				Italian	Non sterile
				Japanese	未滅菌
				Latvian	Nesterils
				Lithuanian	Nesterilus
				Norwegian	Usteril
				Polish	Niesterylne
				Portuguese	Não esterilizado
Romanian	Nesteril				
Slovak	Nesterilné				
Slovenian	Nesterilno				
Spanish	No estéril				
Swedish	Icke-steril				
Turkish	Steril değil				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.4.3	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the need for the user to consult the <i>instructions for use</i>	English	Consult instructions for use
				Bulgarian	Вижте инструкциите за употреба
				Croatian	Pogledajte upute za upotrebu
				Chinese	参阅使用说明书
				Czech	Přečtěte si návod k použití
				Danish	Se brugervejledning
				Dutch	Zie gebruiksaanwijzing
				Estonian	Lugege kasutusjuhendit
				Finnish	Katso käyttöohje
				French	Consulter le mode d'emploi
				German	Gebrauchsanweisung befolgen
				Greek	Συμβουλευτείτε τις οδηγίες χρήσης
				Hungarian	Olvassa el a használati útmutatót
				Italian	Consultare le istruzioni per l'uso
				Japanese	使用方法についてはお問い合わせください
				Latvian	Skatīt lietošanas instrukciju
				Lithuanian	Informacijos ieškokite naudojimo instrukcijoje
				Norwegian	Se bruksanvisningen
				Polish	Należy zapoznać się z instrukcją użytkowania
				Portuguese	Consulte as instruções de utilização
Romanian	Consultați instrucțiunile de utilizare				
Slovak	Prečítajte si návod na použitie				
Slovenian	Preberite navodila za uporabo				
Spanish	Consultar las instrucciones de uso				
Swedish	Se bruksanvisningen				
Turkish	Kullanma kılavuzuna başvurun				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.4.4	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action to avoid undesirable consequences	English	Caution
				Bulgarian	Внимание
				Croatian	Oprez
				Chinese	注意
				Czech	Upozornění
				Danish	Forsigtig
				Dutch	Let op
				Estonian	Ettevaatust!
				Finnish	Varoitus
				French	Attention
				German	Achtung
				Greek	Προσοχή
				Hungarian	Vigyázat!
				Italian	Attenzione
				Japanese	注意
				Latvian	Uzmanību!
				Lithuanian	Perspėjimas
				Norwegian	Forsiktig
				Polish	Ostrożnie
				Portuguese	Cuidado
Romanian	Atenție				
Slovak	Upozornenie				
Slovenian	Pozor				
Spanish	Precaución				
Swedish	Försiktighet				
Turkish	Dikkat				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.1.8	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the entity importing the <i>medical device</i> into a market	English	Importer
				Bulgarian	Вносител
				Croatian	Uvoznik.
				Chinese	进口商
				Czech	Dovozce
				Danish	Importør
				Dutch	Importeur
				Estonian	Importija
				Finnish	Maahantuojaja
				French	Importateur
				German	Importeur
				Greek	Εισαγωγέας
				Hungarian	Importőr
				Italian	Importatore
				Japanese	輸入業者
				Latvian	Importētājs
				Lithuanian	Importuotojas
				Norwegian	Importør
				Polish	Importer
				Portuguese	Importador
Romanian	Importator				
Slovak	Dovozca				
Slovenian	Uvoznik				
Spanish	Importador				
Swedish	Importör				
Turkish	İthalatçı				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.7.10	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements.	Indicates a carrier that contains unique device identifier information.	English	Unique device Identifier
				Bulgarian	Уникален идентификатор на изделията
				Croatian	Jedinstveni identifikator proizvoda
				Chinese	唯一设备标识符
				Czech	Jedinečný identifikátor prostředku
				Danish	Unik udstyrsidentifikation
				Dutch	Unieke hulpmiddelenidentificatie
				Estonian	Kordumatu identifitseerimistunnus
				Finnish	Yksilöllinen laitetunniste
				French	Identifiant unique du dispositif
				German	Eindeutige Produktkennung
				Greek	Μοναδικό αναγνωριστικό συσκευής
				Hungarian	Egyedi eszközazonosító
				Italian	Identificativo univoco del dispositivo
				Japanese	固有のデバイス識別子
				Latvian	Unikālais ierīces identifikators
				Lithuanian	Unikalusis prietaiso identifikatorius
				Norwegian	Unik enhetsidentifikator
				Polish	Niepowtarzalny kod identyfikacyjny wyrobu
				Portuguese	Identificador único do dispositivo
				Romanian	Identificatori unici ai dispozitivelor
				Slovak	Unikátny identifikátor pomôcky
				Slovenian	Edinstveni identifikator pripomočka
Spanish	Identificador único del producto				
Swedish	Unik produktidentifierare				
Turkish	Benzersiz cihaz tanımlayıcısı				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements Combined with the standard form for negation stated in Annex B	Indicates that the <i>medical device</i> is not made with natural rubber latex	English	Not made with natural rubber latex
				Bulgarian	Не е произведено с естествен каучуков латекс
				Croatian	Nije proizvedeno s prirodnim gumenim lateksom
				Chinese	不含天然乳胶
				Czech	Vyrobeno bez použití přírodního kaučukového latexu
				Danish	Ikke fremstillet af naturlig gummilatex
				Dutch	Niet gemaakt met natuurlijke rubber latex
				Estonian	Ei sisalda looduslikku kummilateksit
				Finnish	Valmistuksessa ei ole käytetty luonnonkumilateksia
				French	Ne contient aucun latex de caoutchouc naturel
				German	Ohne Naturkautschuklatex hergestellt
				Greek	Δεν είναι κατασκευασμένο με λάτεξ από φυσικό καουτσούκ
				Hungarian	Nem tartalmaz természetes gumilatekset
				Italian	Non contiene lattice di gomma naturale
				Japanese	天然ゴムラテックスは使用されていません
				Latvian	Nav izgatavots no dabiskā kaučuka lateksa
				Lithuanian	Sudėtyje nėra natūralaus kaučiuko latekso
				Norwegian	Ikke fremstilt med naturgummilateks
				Polish	Nie został wykonany z lateksu z kauczuku naturalnego
				Portuguese	Isento de látex de borracha natural
Romanian	Nu este fabricat cu latex din cauciuc natural				
Slovak	Nie je vyrobený s použitím prírodného gumového latexu				
Slovenian	Ni izdelano iz naravnega lateksa				
Spanish	<i>Fabricado sin látex de caucho natural</i>				
Swedish	Innehåller inte naturgummilatex				
Turkish	Doğal kauçuk lateksten üretilmemiştir				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.4.10	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine disrupting properties	English	Contains hazardous substances
				Bulgarian	Съдържа опасни вещества
				Croatian	Sadrži opasne tvari
				Chinese	含有危险物质
				Czech	Obsahuje nebezpečné látky
				Danish	Indeholder farlige stoffer
				Dutch	Bevat gevaarlijke stoffen
				Estonian	Sisaldab ohtlikke aineid
				Finnish	Sisältää vaarallisia aineita
				French	Contient des substances dangereuses
				German	Enthält gefährliche Stoffe
				Greek	Περιέχει επικίνδυνες ουσίες
				Hungarian	Veszélyes anyagokat tartalmaz
				Italian	Contiene sostanze pericolose
				Japanese	有害物質を含む
				Latvian	Satur bīstamas vielas
				Lithuanian	Sudėtyje yra pavojingų medžiagų
				Norwegian	Inneholder farlige stoffer
				Polish	Zawiera substancje niebezpieczne
				Portuguese	Contém substâncias perigosas
Romanian	Conține substanțe periculoase				
Slovak	Obsahuje nebezpečné látky				
Slovenian	Vsebuje nevarne snovi				
Spanish	Contiene sustancias peligrosas				
Swedish	Innehåller farliga ämnen				
Turkish	Tehlikeli maddeler içerir				



Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.2.8	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates that a <i>medical device</i> should not be used if the package has been damaged or opened, and that the user should consult the <i>instructions for use</i> for additional information	English	Do not use if package is damaged
				Bulgarian	Не употребявайте, ако опаковката е повредена
				Croatian	Nemojte upotrebljavati ako je pakiranje oštećeno
				Chinese	包装损坏时请勿使用
				Czech	Nepoužívejte, pokud je obal poškozen
				Danish	Må ikke anvendes, hvis pakningen er beskadiget
				Dutch	Niet gebruiken als de verpakking beschadigd is
				Estonian	Mitte kasutada, kui pakend on kahjustatud
				Finnish	Käyttö kielletty, jos pakkaus on vaurioitunut
				French	Ne pas utiliser le produit si l'emballage est endommagé
				German	Nicht verwenden, wenn die Verpackung beschädigt ist
				Greek	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά
				Hungarian	Ne használja, ha a csomagolás sérült
				Italian	Non utilizzare se la confezione è danneggiata
				Japanese	梱包が損傷している場合、使用しないでください
				Latvian	Nelietot, ja iepakojums ir bojāts
				Lithuanian	Nenaudoti, jei pakuotė pažeista
				Norwegian	Må ikke brukes hvis pakken er skadet
				Polish	Nie używać, jeśli opakowanie jest uszkodzone
				Portuguese	Não utilize se a embalagem estiver danificada
				Romanian	A nu se utiliza dacă ambalajul este deteriorat
Slovak	Nepoužívajte, ak je balenie poškodené				
Slovenian	Ne uporabljajte, če je embalaža poškodovana				
Spanish	No utilice el dispositivo si el embalaje está dañado				
Swedish	Använd inte utrustningen om förpackningen är skadad.				
Turkish	Paket hasarlıysa kullanmayın				



Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.3.2	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a <i>medical device</i> that needs protection from light sources	English	Keep away from sunlight
				Bulgarian	Да се пази от слънчева светлина
				Croatian	Držite dalje od sunčeve svjetlosti
				Chinese	远离阳光
				Czech	Chraňte před slunečním světlem
				Danish	Må ikke udsættes for sollys
				Dutch	Uit de buurt van zonlicht houden
				Estonian	Kaitsta otsese päikesevalguse eest
				Finnish	Säilytä auringonvalolta suojattuna
				French	Conserver à l'abri de la lumière du soleil
				German	Vor Sonnenlicht schützen
				Greek	Διατηρείτε μακριά από το ηλιακό φως
				Hungarian	Napfénytől elzárva tartandó
				Italian	Tenere lontano dalla luce solare
				Japanese	直射日光を避けてください
				Latvian	Sargāt no saules stariem
				Lithuanian	Saugoti nuo saulės šviesos
				Norwegian	Holdes unna sollys
				Polish	Nie wystawiać na światło słoneczne
				Portuguese	Manter afastado da luz solar
Romanian	A se feri de lumina solară				
Slovak	Chránite pred slnečným žiarením				
Slovenian	Zaščitite pred sončnimi žarki				
Spanish	Mantener alejado de la luz solar				
Swedish	Utsätt inte för direkt sol				
Turkish	Güneş ışığından uzak tutun				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.3.4	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates a medical device that needs to be protected from moisture	English	Keep dry
				Bulgarian	Съхранявайте на сухо
				Croatian	Čuvajte na suhom
				Chinese	保持干燥
				Czech	Uchovávejte v suchu
				Danish	Opbevares tørt
				Dutch	Droog houden
				Estonian	Hoida kuivas
				Finnish	Säilytä kuivassa
				French	Conserver au sec
				German	Trocken aufbewahren
				Greek	Διατηρείτε στεγνό
				Hungarian	Szárazon tartandó
				Italian	Mantenere asciutto
				Japanese	湿気厳禁
				Latvian	Glabāt sausumā
				Lithuanian	Laikyti sausoje vietoje
				Norwegian	Oppbevares tørt
				Polish	Przechowywać w suchym miejscu
				Portuguese	Manter seco
Romanian	A se păstra uscat				
Slovak	Uchovávať v suchu				
Slovenian	Hranite na suhem				
Spanish	Mantener seco				
Swedish	Förvaras torrt				
Turkish	Kuru tutun				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.1.7	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the manufacturer's serial number so that a specific medical device can be identified	English	Serial number
				Bulgarian	Сериен номер
				Croatian	Serijski broj
				Chinese	序列号
				Czech	Sériové číslo
				Danish	Serienummer
				Dutch	Serienummer
				Estonian	Seerianumber
				Finnish	Sarjanumero
				French	Numéro de série
				German	Seriennummer
				Greek	Σειριακός αριθμός
				Hungarian	Sorozatszám
				Italian	Numero di serie
				Japanese	シリアル番号
				Latvian	Sērijas numurs
				Lithuanian	Serijos numeris
				Norwegian	Serienummer
				Polish	Numer seryjny
				Portuguese	Número de série
				Romanian	Număr de serie
Slovak	Sériové číslo				
Slovenian	Serijska številka				
Spanish	Número de serie				
Swedish	Serienummer				
Turkish	Seri numarası				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.7.7	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the item is a <i>medical device</i>	English	Medical device
				Bulgarian	Медицинско изделие
				Croatian	Medicinski proizvod
				Chinese	医疗器械
				Czech	Zdravotnický prostředek
				Danish	Medicinsk udstyr
				Dutch	Medisch hulpmiddel
				Estonian	Meditsiiniseade
				Finnish	Lääkinnällinen laite
				French	Dispositif médical
				German	Medizinprodukt
				Greek	Ιατροτεχνολογικό προϊόν
				Hungarian	Orvostechikai eszköz
				Italian	Dispositivo medico
				Japanese	医療装置
				Latvian	Medicīniskā ierīce
				Lithuanian	Medicinos prietaisas
				Norwegian	Medisinsk enhet
				Polish	Wyrób medyczny
				Portuguese	Dispositivo médico
Romanian	Dispozitiv medical				
Slovak	Zdravotnícka pomôcka				
Slovenian	Medicinski pripomoček				
Spanish	Producto sanitario				
Swedish	Medicinteknisk produkt				
Turkish	Tıbbi cihaz				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 <p>Example:</p> 	5.3.9	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the range of atmospheric pressure to which the <i>medical device</i> can be safely exposed	English	Atmospheric pressure limitation
				Bulgarian	Ограничение на атмосферното налягане
				Croatian	Granične vrijednosti atmosferskog tlaka
				Chinese	大气压强限制
				Czech	Omezení atmosférickým tlakem
				Danish	Atmosfærisk trykgrænse
				Dutch	Bereik van de atmosferische druk
				Estonian	Õhurõhupiirang
				Finnish	Ilmanpaineen raja-arvot
				French	Limites de pression atmosphérique
				German	Beschränkung des atmosphärischen Drucks
				Greek	Περιορισμός ατμοσφαιρικής πίεσης
				Hungarian	Légköri nyomás korlátozása
				Italian	Limite di pressione atmosferica
				Japanese	大気圧制限
				Latvian	Atmosfēras spiediena ierobežojums
				Lithuanian	Atmosferinio slėgio ribos
				Norwegian	Begrensning i atmosfærisk trykk
				Polish	Ograniczenie ciśnienia atmosferycznego
				Portuguese	Limite de pressão atmosférica
Romanian	Limite de presiune atmosferică				
Slovak	Obmedzenie atmosférického tlaku				
Slovenian	Omejitev atmosferskega tlaka				
Spanish	Límites de presión atmosférica				
Swedish	Atmosfärstryckbegränsning				
Turkish	Atmosferik basınç sınırlaması				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 <p>Example:</p> 	5.3.8	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	Indicates the range of humidity to which the <i>medical device</i> can be safely exposed	English	Humidity limitation
				Bulgarian	Ограничение на влажността
				Croatian	Granične vrijednosti vlage
				Chinese	湿度限制
				Czech	Vlhkostní omezení
				Danish	Luffugtighedsgrænse
				Dutch	Vochtigheidsbereik
				Estonian	Niiskuse piirang
				Finnish	Kosteuden raja-arvot
				French	Limites d'humidité
				German	Feuchtigkeitsbeschränkung
				Greek	Περιορισμός υγρασίας
				Hungarian	Páratartalom-korlátozás
				Italian	Limite di umidità
				Japanese	湿度制限
				Latvian	Mitruma ierobežojums
				Lithuanian	Drėgmės ribos
				Norwegian	Fuktighetsbegrensning
				Polish	Ograniczenie wilgotności
				Portuguese	Limite de humidade
Romanian	Limite de umiditate				
Slovak	Obmedzenie vlhkosti				
Slovenian	Omejitev vlažnosti				
Spanish	Límites de humedad				
Swedish	Fuktighetsbegränsning				
Turkish	Nemlilik sınırlaması				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5.1.2	ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements	<p>This symbol is replaced by the XX-REP symbol as per ISO 15223-1:2021. The EC-REP symbol is kept in this document until all devices with this symbol are no longer in use.</p> <p>Indicates the authorized representative in the European Community/European Union</p>	English	Authorised representative in the European Community/European Union
				Bulgarian	Упълномощен представител в Европейската общност/Европейския съюз
				Croatian	Ovlašteni predstavnik u Europskoj zajednici/Europskoj uniji
				Chinese	欧洲共同体 / 欧盟授权代表
				Czech	Autorizovaný zástupce v Evropském společenství/Evropské unii
				Danish	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union
				Dutch	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/Europese Unie
				Estonian	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus
				Finnish	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa
				French	Mandataire dans la Communauté européenne/l'Union européenne
				German	Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union
				Greek	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/Ευρωπαϊκή Ένωση
				Hungarian	Meghatalmazott képviselő az Európai Közösségekben/Európai Unióban
				Italian	Rappresentante autorizzato nella Comunità europea/Unione europea
Japanese	欧州共同体 / 欧州連合における正式代理人				


				Latvian	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā
				Lithuanian	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
				Norwegian	Autorisert representant i EU/EØS
				Polish	Autoryzowany przedstawiciel we Wspólnocie Europejskiej/Unii Europejskiej
				Portuguese	Representante autorizado na Comunidade Europeia/União Europeia
				Romanian	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană
				Slovak	Autorizovaný zástupca v Európskom spoločenstve/Európskej únii
				Slovenian	Pooblaščen zastopnik v Evropski skupnosti/Evropski uniji
				Spanish	Representante autorizado en la Comunidad Europea o la Unión Europea
				Swedish	Auktoriserad representant i Europeiska gemenskapen/Europeiska unionen
				Turkish	Avrupa Topluluğu/Avrupa Birliği yetkili temsilcisi

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title		
<table border="1"> <tr> <td>XX</td> <td>REP</td> </tr> </table>	XX	REP	5.1.2	ISO 15223-1:2021/Amd 1:2025 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements AMENDMENT 1: Addition of defined	This symbol replaces the EC-REP symbol, see above. Indicates the <i>authorized representative</i> in the identified country or jurisdiction.	English	Authorized representative in the identified country or jurisdiction
	XX	REP					
	Bulgarian	Оторизиран представител в определената държава или юрисдикция					
	Croatian	Ovlašteni predstavnik u određenoj zemlji ili jurisdikciji					
	Chinese	指定国家或地区的授权代表					
Czech	Zplnomocněný zástupce v příslušné zemi nebo jurisdikci						

		term for authorized representative and modified EC REP symbol to not be country or region specific	The [XX] text of the <i>symbol</i> shall be replaced by either the two-letter country code or the three-letter country code defined in ISO 3166-1 or other text required by authority having jurisdiction. This <i>symbol</i> shall be accompanied by the name and address of the <i>authorized representative</i> adjacent to the <i>symbol</i> .	Danish	Autoriseret repræsentant i det identificerede land eller den identificerede jurisdiktion
				Dutch	Gemachtigde vertegenwoordiger in het geïdentificeerde land of rechtsgebied
				Estonian	Volitatud esindaja määratletud riigis või jurisdiktsioonis
				Finnish	Valtuutettu edustaja määritetyssä maassa tai lainkäyttöalueella
				French	Représentant agréé dans le pays ou la juridiction identifiés
				German	Autorisierte Vertretung in dem jeweiligen Land oder der jeweiligen Gerichtsbarkeit
				Greek	Εξουσιοδοτημένος αντιπρόσωπος στην προσδιορισμένη χώρα ή δικαιοδοσία
				Hungarian	Meghatalmazott képviselő a megjelölt országban vagy joghatóságban
				Italian	Rappresentante autorizzato per il paese o la giurisdizione identificati
				Japanese	指定された国または管轄区域における認定代表者
				Latvian	Pilnvarotais pārstāvis norādītajā valstī vai jurisdikcijā
				Lithuanian	Įgaliotasis atstovas nurodytoje šalyje arba jurisdikcijoje
				Norwegian	Autorisert representant i det identifiserte landet eller jurisdiksjonen
				Polish	Autoryzowany przedstawiciel w określonym kraju lub jurysdykcji
				Portuguese	Representante autorizado no país ou jurisdição identificados
				Romanian	Reprezentant autorizat în țara sau jurisdicția identificată
				Slovak	Spĺnomocnený zástupca v určenej krajine alebo jurisdikcii
				Slovenian	Pooblaščeni predstavnik v določeni državi ali jurisdikciji
				Spanish	Representante autorizado en el país o jurisdicción identificados
				Swedish	Auktoriserad representant i angivet land eller angiven jurisdiktion
				Turkish	Belirtilen ülke veya yargı bölgesindeki yetkili temsilci


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	ISO 15986 Symbols for use in the labelling of medical devices – Requirements for labelling of medical devices containing phthalates Combined with the standard form for negation stated in Annex B	Indicates that the <i>medical device</i> does not contain phthalates	English	Product is phthalates free
				Bulgarian	Продуктът не съдържа фталати
				Croatian	Proizvod ne sadrži ftalate
				Chinese	产品不含邻苯二甲酸盐
				Czech	Výrobek neobsahuje ftaláty
				Danish	Produktet er phthalatfrit
				Dutch	Product is ftalaatvrij
				Estonian	Toode ei sisalda ftalaate
				Finnish	Tuote ei sisällä ftalaatteja
				French	Absence de phtalate
				German	Produkt ist frei von Phthalaten
				Greek	Το προϊόν δεν περιέχει φθαλικές ενώσεις
				Hungarian	A termék ftalátmentes
				Italian	Il prodotto non contiene ftalati
				Japanese	製品にはフタル酸エステル類は使用されていません
				Latvian	Izstrādājums nesatur ftalātus
				Lithuanian	Gaminyje nėra ftalatų
				Norwegian	Produktet inneholder ikke ftalater
				Polish	Produkt nie zawiera ftalanów
				Portuguese	O produto não contém ftalatos
Romanian	Produsul nu conține ftalați				
Slovak	Výrobok neobsahuje ftaláty				
Slovenian	Izdelek ne vsebuje ftalotov				
Spanish	El producto no contiene ftalatos				
Swedish	Produkten är ftalatfri				
Turkish	Ürün ftalat içermez				


MDR 2017/745


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	Medical Device Regulation 2017/745	Indicates that the product is in compliance with European legislation for medical devices	English	European Conformity
				Bulgarian	Европейско съответствие
				Croatian	Sukladnost s europskim normama
				Chinese	欧盟认证
				Czech	Shoda s právními předpisy Evropské unie
				Danish	Europæisk overensstemmelse
				Dutch	Europese conformiteit
				Estonian	Euroopa vastavus
				Finnish	CE-merkintä, Euroopan lainsäädännön mukainen
				French	Marquage CE
				German	Europäische Konformität
				Greek	Ευρωπαϊκή συμμόρφωση
				Hungarian	Európai megfeleléség
				Italian	Conformità europea
				Japanese	欧州適合性
				Latvian	Eiropas atbilstības marķējums
				Lithuanian	Europos atitiktis
				Norwegian	Europeisk samsvar
				Polish	Zgodność europejska
				Portuguese	Conformidade europeia
Romanian	Conformitate europeană				
Slovak	Zhoda s európskymi normami				
Slovenian	Evropska skladnost				
Spanish	Conformidad europea				
Swedish	Överensstämmelse med EU-standard				
Turkish	Avrupa Uygunluğu				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	Medical Device Regulation 2017/745	Indicates that the product is in compliance with European legislation for medical devices	English	European Conformity
				Bulgarian	Европейско съответствие
				Croatian	Sukladnost s europskim normama
				Chinese	欧盟认证
				Czech	Shoda s právními předpisy Evropské unie
				Danish	Europæisk overensstemmelse
				Dutch	Europese conformiteit
				Estonian	Euroopa vastavus
				Finnish	CE-merkintä, Euroopan lainsäädännön mukainen
				French	Marquage CE
				German	Europäische Konformität
				Greek	Ευρωπαϊκή συμμόρφωση
				Hungarian	Európai megfeleléség
				Italian	Conformità europea
				Japanese	欧州適合性
				Latvian	Eiropas atbilstības marķējums
				Lithuanian	Europos atitiktis
				Norwegian	Europeisk samsvar
				Polish	Zgodność europejska
				Portuguese	Conformidade europeia
				Romanian	Conformitate europeană
Slovak	Zhoda s európskymi normami				
Slovenian	Evropska skladnost				
Spanish	Conformidad europea				
Swedish	Överensstämmelse med EU-standard				
Turkish	Avrupa Uygunluğu				


Medicines and Medical Devices Act 2021 – UK only


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	The Medical Devices Regulations 2002 (SI 2002 No 618) as subsequently amended by the EU Exit Regulations of 2019 (SI 791) and 2020 (SI 1478).	Indicates that a product is compliant with British legislation for medical devices	English	UK Conformity Assessed
				Bulgarian	Оценено за съответствие за Обединеното кралство
				Croatian	Procjena sukladnosti za UK
				Chinese	英国合格认定
				Czech	Posouzení shody s předpisy Velké Británie
				Danish	Storbritanniens overensstemmelse vurderet
				Dutch	VK-conformiteit beoordeeld
				Estonian	Ühendkuningriigi vastavushinnanguga
				Finnish	UKCA-merkintä, Yhdistyneen kuningaskunnan lainsäädännön mukainen
				French	Conformité UKCA
				German	UK-Konformitätsbewertung
				Greek	Με πιστοποίηση συμμόρφωσης του ΗΒ
				Hungarian	Egyesült királysági megfelelésértékelve
				Italian	Conformità Regno Unito valutata
				Japanese	英国適合性評価済み
				Latvian	AK atbilstība novērtēta
				Lithuanian	JK atitiktis įvertinta
				Norwegian	UK-samsvar vurdert
				Polish	Ocena zgodności w Wielkiej Brytanii
				Portuguese	Avaliação da conformidade do Reino Unido
Romanian	Conformitate UK evaluată				
Slovak	Vyhodnotená zhoda s normami Spojeného kráľovstva				
Slovenian	Ocenjena skladnost ZK				
Spanish	Evaluación británica de la conformidad				
Swedish	Bedömning av överensstämmelse i Storbritannien				
Turkish	Birleşik Krallık Uygunluk Değerlendirmesi				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	The Medical Devices Regulations 2002 (SI 2002 No 618) as subsequently amended by the EU Exit Regulations of 2019 (SI 791) and 2020 (SI 1478).	Indicates that a product is compliant with British legislation for medical devices	English	UK Conformity Assessed
				Bulgarian	Оценено за съответствие за Обединеното кралство
				Croatian	Procjena sukladnosti za UK
				Chinese	英国合格认定
				Czech	Posouzení shody s předpisy Velké Británie
				Danish	Storbritanniens overensstemmelse vurderet
				Dutch	VK-conformiteit beoordeeld
				Estonian	Ühendkuningriigi vastavushinnangu
				Finnish	UKCA-merkintä, Yhdistyneen kuningaskunnan lainsäädännön mukainen
				French	Conformité UKCA
				German	UK-Konformitätsbewertung
				Greek	Με πιστοποίηση συμμόρφωσης του ΗΒ
				Hungarian	Egyesült királysági megfelelésértékelve
				Italian	Conformità Regno Unito valutata
				Japanese	英国適合性評価済み
				Latvian	AK atbilstība novērtēta
				Lithuanian	JK atitiktis įvertinta
				Norwegian	UK-samsvar vurdert
				Polish	Ocena zgodności w Wielkiej Brytanii
				Portuguese	Avaliação da conformidade do Reino Unido
Romanian	Conformitate UK evaluată				
Slovak	Vyhodnotená zhoda s normami Spojeného kráľovstva				
Slovenian	Ocenjena skladnost ZK				
Spanish	Evaluación británica de la conformidad				
Swedish	Bedömning av överensstämmelse i Storbritannien				
Turkish	Birleşik Krallık Uygunluk Değerlendirmesi				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	The Medical Devices Regulations 2002 (SI 2002 No 618) as subsequently amended by the EU Exit Regulations of 2019 (SI 791) and 2020 (SI 1478).	Indicates an UK Responsible person	English	UK Responsible Person
				Bulgarian	Отговорно лице за Обединеното кралство
				Croatian	Odgovorna osoba za Ujedinjenu Kraljevinu
				Chinese	英国负责人
				Czech	Odpovědná osoba ve Spojeném království
				Danish	Ansvarshavende i UK
				Dutch	Verantwoordelijke voor het VK
				Estonian	Vastutav isik Ühendkuningriigis
				Finnish	Vastuhenkilö Yhdistyneessä kuningaskunnassa
				French	Responsable Royaume-Uni
				German	Verantwortliche Person im Vereinigten Königreich
				Greek	Αρμόδιο πρόσωπο στο Ηνωμένο Βασίλειο
				Hungarian	Egyesült királysági felelős személy
				Italian	Persona responsabile nel Regno Unito
				Japanese	英国責任者
				Latvian	Atbildīgā persona AK
				Lithuanian	JK atsakingas asmuo
				Norwegian	Ansvarlig person i Storbritannia
				Polish	Osoba odpowiedzialna w Wielkiej Brytanii
				Portuguese	Pessoa responsável no Reino Unido
				Romanian	Persoana responsabilă în Marea Britanie
Slovak	Zodpovedná osoba pre Spojené kráľovstvo				
Slovenian	Odgovorna oseba v Združenem kraljestvu				
Spanish	Persona responsable en el Reino Unido				
Swedish	Ansvarig person, Storbritannien				
Turkish	BK Sorumlusu				

Directive 2012/19/EU – EN 50419 (E)

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE) EN 50419:2022 (E) – Marking of electrical and electronic equipment (EEE) in respect to separate of waste EEE (WEEE)	Indicates that separate collection for electrical and electronic equipment is required according to local regulation and collection schemes for disposal of electronic and electrical waste (WEEE)	English	Separate collection of electrical and electronic equipment
				Bulgarian	Разделно събиране на електрическо и електронно оборудване
				Croatian	Odvojeno prikupljanje električne i elektroničke opreme
				Chinese	单独收集电气和电子设备
				Czech	Oddělený sběr elektrických a elektronických zařízení
				Danish	Separat indsamling af elektrisk og elektronisk udstyr
				Dutch	Gescheiden inzameling van elektrische en elektronische apparatuur
				Estonian	Elektri- ja elektroonikaseadmete eraldi kogumine
				Finnish	Sähkö- ja elektroniikkalaiteromun erilliskeräys
				French	Collecte séparée des équipements électriques et électroniques
				German	Getrennte Sammlung von Elektro- und Elektronikgeräten
				Greek	Ξεχωριστή συλλογή ηλεκτρικού και ηλεκτρονικού εξοπλισμού
				Hungarian	Az elektromos és elektronikus berendezések hulladékainak szelektív gyűjtése
				Italian	Raccolta differenziata di apparecchiature elettriche ed elettroniche
				Japanese	電気・電子機器の分別回収
				Latvian	Atsevišķa elektrisko un elektronisko iekārtu savākšana
				Lithuanian	Atskiras elektros ir elektroninės įrangos surinkimas
				Norwegian	Separat innsamling av elektrisk og elektronisk utstyr
				Polish	Oddzielna selektywna zbiórka sprzętu elektrycznego i elektronicznego
				Portuguese	Recolha diferenciada de equipamento elétrico e eletrónico
				Romanian	Colectarea separată a echipamentelor electrice și electronice
				Slovak	Separovaný zber elektrických a elektronických zariadení
				Slovenian	Ločeno zbiranje električne in elektronske opreme
				Spanish	Recogida selectiva de aparatos eléctricos y electrónicos
Swedish	Separat insamling av elektrisk och elektronisk utrustning				
Turkish	Elektrikli ve elektronik ekipmanların ayrı toplanması				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE) EN 50419:2022 (E) – Marking of electrical and electronic equipment (EEE) in respect to separate collection of waste EEE (WEEE)	Indicates that separate collection for electrical and electronic equipment is required according to local regulation and collection schemes for disposal of electronic and electrical waste (WEEE)	English	Separate collection of electrical and electronic equipment
				Bulgarian	Разделно събиране на електрическо и електронно оборудване
				Croatian	Odvojeno prikupljanje električne i elektroničke opreme
				Chinese	单独收集电气和电子设备
				Czech	Oddělený sběr elektrických a elektronických zařízení
				Danish	Separat indsamlng af elektrisk og elektronisk udstyr
				Dutch	Gescheiden inzameling van elektrische en elektronische apparatuur
				Estonian	Elektri- ja elektroonikaseadmete eraldi kogumine
				Finnish	Sähkö- ja elektroniikkalaiteromun erilliskeräys
				French	Collecte séparée des équipements électriques et électroniques
				German	Getrennte Sammlung von Elektro- und Elektronikgeräten
				Greek	Ξεχωριστή συλλογή ηλεκτρικού και ηλεκτρονικού εξοπλισμού
				Hungarian	Az elektromos és elektronikus berendezések hulladékainak szelektív gyűjtése
				Italian	Raccolta differenziata di apparecchiature elettriche ed elettroniche
				Japanese	電気・電子機器の分別回収
				Latvian	Atsevišķa elektrisko un elektronisko iekārtu savākšana
				Lithuanian	Atskiras elektros ir elektroninės įrangos surinkimas
				Norwegian	Separat innsamlng av elektrisk og elektronisk utstyr
				Polish	Oddzielna selektywna zbiórka sprzętu elektrycznego i elektronicznego
				Portuguese	Recolha diferenciada de equipamento eléctrico e electrónico
				Romanian	Colectarea separată a echipamentelor electrice și electronice
Slovak	Separovaný zber elektrických a elektronických zariadení				
Slovenian	Ločeno zbiranje električne in elektronske opreme				
Spanish	Recogida selectiva de aparatos eléctricos y electrónicos				
Swedish	Separat insamlng av elektrisk och elektronisk utrustning				
Turkish	Elektrikli ve elektronik ekipmanların ayrı toplanması				

EN 60529+A1+A2+AC 2019


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
IP30	N/A	EN 60529+A1 - Degrees of protection provided by enclosures (IP Code)	Degrees of protection provided by enclosures (IP Code)	English	Protection against solid objects
				Bulgarian	Защита срещу твърди предмети
				Croatian	Zaštita od čvrstih predmeta
				Chinese	防止固体碰撞
				Czech	Ochrana proti vniknutí pevných těles
				Danish	Beskyttelse mod faste genstande
				Dutch	Bescherming tegen massieve voorwerpen
				Estonian	Kaitse tahkete esemete eest
				Finnish	Suojaus kiinteitä esineitä vastaan
				French	Protection contre les objets solides
				German	Schutz gegen feste Fremdkörper
				Greek	Προστασία έναντι συμπαγών αντικειμένων
				Hungarian	Védelem szilárd tárgyak ellen
				Italian	Protezione contro gli oggetti solidi
				Japanese	ソリッドオブジェクトへの耐性仕様
				Latvian	Aizsardzība pret cietiem priekšmetiem
				Lithuanian	Apsauga nuo kietų objektų
				Norwegian	Beskyttelse mot faste gjenstander
				Polish	Ochrona przed ciałami stałymi
				Portuguese	Proteção contra objetos sólidos
Romanian	Protecție împotriva obiectelor solide				
Slovak	Ochrana proti pevným časticiam				
Slovenian	Zaščita pred trdnimi predmeti				
Spanish	Protección contra objetos sólidos				
Swedish	Skydd mot fasta föremål				
Turkish	Katı nesnelere karşı koruma				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
IP31	N/A	EN 60529+A1 - Degrees of protection provided by enclosures (IP Code)	Degrees of protection provided by enclosures (IP Code)	English	Protection against solid particles and liquid ingress.
				Bulgarian	Защита срещу твърди частици и проникване на течности.
				Croatian	Zaštita od ulaska krutih čestica i tekućine.
				Chinese	防止固体颗粒和液体进入。
				Czech	Ochrana před vniknutím pevných cizích těles a kapalin.
				Danish	Beskyttelse mod faste partikler og væskeindtrængen.
				Dutch	Bescherming tegen het binnendringen van vaste deeltjes en vloeistoffen.
				Estonian	Kaitse tahkete osakeste ja vedeliku sissetungimise eest
				Finnish	Suojaus kiinteitä hiukkasia ja nesteen tunkeutumista vastaan.
				French	Protection contre les particules solides et la pénétration de liquides.
				German	Schutz gegen Feststoffpartikel und Eindringen von Flüssigkeiten.
				Greek	Προστασία από στερεά σωματίδια και διείσδυση υγρών.
				Hungarian	Szilárd részecskékkel és folyadékbehatolással szembeni védelem.
				Italian	Protezione dall'ingresso di particelle solide e liquidi.
				Japanese	固体粒子および液体の侵入に対する保護。
				Latvian	Aizsardzība pret cietām daļiņām un šķidruma iekļūšanu korpusā.
				Lithuanian	Apsauga nuo kietųjų dalelių ir skysčių įsiskverbimo.
				Norwegian	Beskyttelse mot inntrengning av faste partikler og væske.
				Polish	Ochrona przed cząstkami stałymi i wnikaniem cieczy.
				Portuguese	Proteção contra partículas sólidas e infiltração de líquidos.
Romanian	Protecție împotriva particulelor solide și a infiltrării lichidelor.				
Slovak	Ochrana pred pevnými časticami a vniknutím tekutín.				
Slovenian	Zaščitna pred trdnimi delci in vdorom tekočine.				
Spanish	Protección contra la entrada de partículas sólidas y líquidos.				
Swedish	Skydd mot inträngande av fasta föremål och vätskedroppar.				
Turkish	Katı parçacıklara ve sıvı girişine karşı koruma.				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
IP34	N/A	EN 60529+A1 - Degrees of protection provided by enclosures (IP Code)	Protection against solid objects larger than 2.5 mm and against water splashing from any direction.	English	Protection against solid foreign objects and water ingress
				Bulgarian	Защита срещу проникване на твърди чужди тела и вода.
				Croatian	Zaštita od ulaska krutih stranih tijela i vode.
				Chinese	防止固体异物和水进入。
				Czech	Ochrana proti vniknutí pevných cizích předmětů a vniknutí vody.
				Danish	Beskyttelse mod faste fremmedlegemer og væskeindtrængen.
				Dutch	Bescherming tegen vaste vreemde voorwerpen en indringing van water.
				Estonian	Kaitse tahkete võõrkehade ja vee sissetungimise eest.
				Finnish	Suojaus kiinteitä vierasesineitä ja veden tunkeutumista vastaan.
				French	Protection contre la pénétration de corps étrangers solides et d'eau.
				German	Schutz gegen feste Fremdkörper und Eindringen von Wasser.
				Greek	Προστασία από στερεά ξένα σώματα και εισχώρηση νερού.
				Hungarian	Szilárd idegen tárgyak és víz behatolásával szembeni védettség.
				Italian	Protezione contro la penetrazione di corpi estranei solidi e acqua.
				Japanese	固体異物および水の侵入に対する保護。
				Latvian	Aizsardzība pret cietu svešķermeņu un ūdens iekļūšanu korpusā.
				Lithuanian	Apsauga nuo kietųjų svetimkūnių ir vandens įsiskverbimo.
				Norwegian	Beskyttelse mot inntrengning av fremmedlegemer og væske.
				Polish	Ochrona przed wnikaniem ciał stałych i wody.
				Portuguese	Proteção contra corpos estranhos sólidos e entrada de água.
Romanian	Protecție împotriva pătrunderii corpurilor străine solide și a apei.				
Slovak	Ochrana pred vniknutím pevných cudzích predmetov a vniknutím vody.				
Slovenian	Zaščita pred vdorom trdnih tujkov in vode.				
Spanish	Protección contra la penetración de cuerpos sólidos extraños y agua.				
Swedish	Skydd mot inträngande av fasta föremål och vatten.				
Turkish	Katı yabancı cisimler ve su girişine karşı koruma.				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
IP50	N/A	EN 60529+A1 - Degrees of protection provided by enclosures (IP Code)	Degrees of protection provided by enclosures (IP Code)	English	Protection against dust
				Bulgarian	Защита срещу прах
				Croatian	Zaštita od prašine
				Chinese	防尘
				Czech	Ochrana proti prachu
				Danish	Støvbeskyttelse
				Dutch	Bescherming tegen stof
				Estonian	Kaitse tolmu eest
				Finnish	Suojaus pölyä vastaan
				French	Protection contre la poussière
				German	Schutz gegen Staub
				Greek	Προστασία έναντι σκόνης
				Hungarian	Védelem por ellen
				Italian	Protezione contro la polvere
				Japanese	防塵仕様
				Latvian	Aizsardzība pret putekļiem
				Lithuanian	Apsauga nuo dulkių
				Norwegian	Beskyttelse mot støv
				Polish	Ochrona przed pyłami
				Portuguese	Proteção contra poeiras
Romanian	Protecție împotriva prafului				
Slovak	Ochrana proti prachu				
Slovenian	Protiprašna zaščita				
Spanish	Protección contra polvo				
Swedish	Dammskydd				
Turkish	Toza karşı koruma				


ISO 7000/IEC 60417


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5333	ISO 7000:2019 - Graphical symbols for use on equipment – Registered symbols ISO 7000 / IEC Graphical Symbols for use on equipment	Identifies a type BF-applied part	English	Electrical safety type BF applied part
				Bulgarian	Контактна част тип BF за електрическа безопасност
				Croatian	Primijenjen dio električne sigurnosti tipa BF
				Chinese	适用电气安全类型 BF 型应用部分
				Czech	Elektrická bezpečnost: příložná část typu BF
				Danish	Elektrisk sikkerhed, type BF, patientdel
				Dutch	Elektrische veiligheid type BF toegepast onderdeel
				Estonian	Elektriohutus, BF-tüüpi kontaktosa
				Finnish	Sähköturvallisuus: tyyppin BF liityntäosa
				French	Protection contre les chocs électriques: partie appliquée de type BF
				German	Elektrische Sicherheit – Anwendungsteil vom Typ BF
				Greek	Εξάρτημα στο οποίο εφαρμόζεται ηλεκτρική ασφάλεια τύπου BF
				Hungarian	Elektromos bíztonság – BF típusú beteggel érintkező alkatrész
				Italian	Grado di protezione contro shock elettrici: apparecchio di tipo BF
				Japanese	電気安全 BF 型装着部
				Latvian	Elektrodrošības BF tipa daļa, ko izmanto tiešā saskarē ar pacientu
				Lithuanian	Elektros saugos BF tipo besiliečiančioji dalis
				Norwegian	Elektrisk sikkerhet, type BF anvendt del
				Polish	Część zgodna ze standardem typu BF dotyczącym bezpieczeństwa elektrycznego
				Portuguese	Peça aplicada do tipo BF de segurança elétrica
Romanian	Componentă aplicată de tip BF pentru siguranță electrică				
Slovak	Elektrická bezpečnosť aplikovanej časti typu BF				
Slovenian	Del z uporabljenjo električno zaščito tipa BF				
Spanish	Seguridad eléctrica: pieza aplicada tipo BF				
Swedish	Elsäkerhet: typ BF, applicerad del				
Turkish	Elektrik güvenliği tipi BF uygulanan parça				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5334	ISO 7000:2019 - Graphical symbols for use on equipment – Registered symbols ISO 7000 / IEC Graphical Symbols for use on equipment	Identifies a type BF-applied part that is defibrillator-proof	English	Defibrillation proof type BF applied part
				Bulgarian	Контактна част тип BF, устойчива на дефибрилация
				Croatian	Dio s primijenjenom zaštitom tipa BF koji je ispitan na defibrilaciju
				Chinese	适用非除颤类型 BF 型应用部分
				Czech	Příložná část typu BF, odolná proti defibrilátoru
				Danish	Defibrilleringssikker type BF anvendt del
				Dutch	Bestand tegen defibrillatie type BF toegepast onderdeel
				Estonian	Defibrilleerimiskindel BF-tüüpi kontaktosa
				Finnish	Defibrillaatiota kestävä tyypin BF liityntäosa
				French	Partie appliquée de type BF protégée contre les chocs de défibrillation
				German	Defibrillationsschutz – Anwendungsteil vom Typ BF
				Greek	Εφαρμοζόμενο εξάρτημα τύπου BF ανθεκτικό στην απινίδωση
				Hungarian	BF típusú, defibrillátor-biztos, beteggel érintkező alkatrész
				Italian	Parte applicata di tipo BF a prova di defibrillazione
				Japanese	耐除細動型 BF 装着部
				Latvian	Pret defibrilāciju droša BF tipa daļa, ko izmanto tiešā saskarē ar pacientu
				Lithuanian	Defibriliacijai atspari BF tipo besiliečiančioji dalis
				Norwegian	Defibrilleringssikker, type BF anvendt del
				Polish	Odporna na defibrylację część aplikacyjna typu BF
				Portuguese	Peça aplicada do tipo BF à prova de desfibrilhação
Romanian	Componentă aplicată de tip BF rezistentă la defibrilare				
Slovak	Aplikovaná časť typu BF odolná voči defibrilácii				
Slovenian	Del z uporabljeno električno zaščito tipa BF, preizkušen na defibrilacijo				
Spanish	Pieza aplicada tipo BF a prueba de desfibrilación				
Swedish	Defibrilleringssäkerhet: typ BF, applicerad del				
Turkish	Defibrilasyona karşı dayanıklılık tipi BF uygulanan parça				


ASTM F2503 – 23


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	ASTM F2503 - 23 Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment	Indicates that the <i>medical device</i> does not pose any safety hazard in MR environments.	English	MR conditional
				Bulgarian	Безопасно за МР при определени условия
				Croatian	Uvjetno sigurno za upotrebu uz magnetsku rezonanciju (MR)
				Chinese	MR 特定条件下安全
				Czech	Podmíněně kompatibilní s MR
				Danish	MR-betinget
				Dutch	MR-voorwaardelijk
				Estonian	MR-tingimuslik
				Finnish	MK-ehdollinen
				French	Compatible IRM sous conditions
				German	Bedingt MR-sicher
				Greek	MR υπό όρους
				Hungarian	MR-környezetben feltételesen biztonságos
				Italian	Compatibilità RM condizionata
				Japanese	MRI には制限付きで対応
				Latvian	Saderīgs ar MR
				Lithuanian	Sąlyginai saugus naudoti MR aplinkoje
				Norwegian	MR-betinget
				Polish	Warunkowo dopuszczzone do stosowania w środowisku rezonansu magnetycznego
				Portuguese	Condicional para RM
Romanian	Compatibilitate condiționată în mediul RM				
Slovak	Podmienečne bezpečný v prostredí MR				
Slovenian	Pogojno varno za MR				
Spanish	Compatibilidad condicional con RM				
Swedish	MR-villkorad				
Turkish	MR koşullu				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	ASTM F2503 - 23 Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment	Indicates that the device poses no known health hazards resulting from exposure to any MR environment	English	MR safe
				Bulgarian	Безопасно за работа в магнитнорезонансна среда
				Croatian	Sigurno za upotrebu uz magnetsku rezonanciju (MR)
				Chinese	MR 安全
				Czech	Bezpečné při použití MR
				Danish	MR-sikker
				Dutch	MR-veilig
				Estonian	MR-ohutu
				Finnish	MK-turvallinen
				French	Compatible IRM
				German	MR-sicher
				Greek	Ασφαλές για MR
				Hungarian	MRI szempontjából biztonságos
				Italian	Compatibile con RM
				Japanese	MR 適合
				Latvian	MR drošs
				Lithuanian	Saugus naudoti MR aplinkoje
				Norwegian	MR-sikker
				Polish	Bezpieczny w trakcie badania rezonansem magnetycznym
				Portuguese	Seguro para RM
				Romanian	Sigur în mediul RM
				Slovak	Bezpečný v prostredí MR
				Slovenian	Varno za MR
Spanish	Compatible con RM				
Swedish	MR-säker				
Turkish	MR açısından güvenli				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	ASTM F2503 - 23 Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment	Indicates that the device poses unacceptable risks to the patient, staff or other persons within the MR environment	English	MR unsafe
				Bulgarian	Не е безопасно за работа в магнитнорезонансна среда
				Croatian	Nije sigurno za upotrebu uz magnetsku rezonanciju (MR)
				Chinese	MR 不安全
				Czech	Nebezpečné při použití MR
				Danish	Ikke MR-sikker
				Dutch	MR-onveilig
				Estonian	Ei ole MR-ohutu
				Finnish	MK-vaarallinen
				French	Non compatible IRM
				German	MR-unsicher
				Greek	Μη ασφαλές για MR
				Hungarian	MRI szempontjából nem biztonságos
				Italian	Non sicuro per RMI
				Japanese	MR 非対応
				Latvian	Nedrošs lietošanai MR vidē
				Lithuanian	Nesaugus naudoti MR aplinkoje
				Norwegian	Ikke MR-sikker
				Polish	Nie jest bezpieczny w środowisku MR
				Portuguese	Inseguro para RM
				Romanian	Nu este sigur în mediul RM
Slovak	Nevhodný do prostredia MR				
Slovenian	Ni varno za MR				
Spanish	Incompatible con RM				
Swedish	Ej MR-säker				
Turkish	MR açısından güvenli değil				



ISO 7010


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	W001	ISO 7010 – Graphical symbols – Safety colours and safety signs EN 60601-1 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.	This sign indicates a hazard with a medium level of risk which, if not avoided, could result in death or serious injury.	English	General warning signs
				Bulgarian	Общи предупредителни знаци
				Croatian	Opći znakovi upozorenja
				Chinese	一般警告标志
				Czech	Obecné výstražné značky
				Danish	Generelle advarsler
				Dutch	Algemene waarschuwborden
				Estonian	Üldised hoiatused
				Finnish	Yleiset varoituserkit
				French	Signaux d'avertissement généraux
				German	Allgemeine Warnzeichen
				Greek	Γενικές σημάνσεις προειδοποίησης
				Hungarian	Általános figyelmeztető jelzések
				Italian	Segnali di avvertimento generale
				Japanese	一般警告フラグ
				Latvian	Vispārējās brīdinājuma zīmes
				Lithuanian	Bendrieji įspėjamieji ženklai
				Norwegian	Generelle advarsler
				Polish	Ogólne znaki ostrzegawcze
				Portuguese	Sinais de aviso gerais
Romanian	Semne de avertizare generală				
Slovak	Všeobecné výstražné symboly				
Slovenian	Splošni opozorilni znaki				
Spanish	Señales de advertencia generales				
Swedish	Allmänna varningssymboler				
Turkish	Genel uyarı işaretleri				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	M002	EN ISO 7010 – Graphical symbols – Safety colours and safety signs – Registered safety signs (ISO 7010:2019) EN 60601-1 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance ISO 20417 Medical devices – Information to be supplied by the manufacturer	Indicates that, according to the manufacturer's risk analysis, not reading the instructions can cause a hazard	English	Action Safety Sign
				Bulgarian	Знак за мерки за безопасност
				Croatian	Sigurnosni znak obveznog postupanja
				Chinese	操作安全标志
				Czech	Čti návod
				Danish	Sikkerhedsskilt
				Dutch	Actieveiligheidssymbool
				Estonian	Tegevuse ohutusmärk
				Finnish	Huomioi turvallisuustoimenpiteet
				French	Symbole de mesure de sécurité
				German	Sicherheitszeichen
				Greek	Σήμανση ασφαλείας σε ενέργεια
				Hungarian	Műveletbiztonsági jel
				Italian	Segnale azione di sicurezza
				Japanese	行動安全表示
				Latvian	Droša rīcība
				Lithuanian	Veiksmų saugos ženklas
				Norwegian	Tiltakssikkerhetsskilt
				Polish	Znaki bezpieczeństwa
				Portuguese	Sinal de segurança da ação
Romanian	Semn pentru măsuri de siguranță				
Slovak	Symbol bezpečnostného opatrenia				
Slovenian	Varnostni znak				
Spanish	Señal de seguridad durante la acción				
Swedish	Säkerhetsskylt				
Turkish	Eylem Güvenlik İşareti				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	M001	EN ISO 7010 – Graphical symbols – Safety colours and safety signs – Registered safety signs IEC 60601-1	Indicates a mandatory action	English	General mandatory action sign
				Bulgarian	Знак за общо задължително действие
				Croatian	Opći znak obveznog postupanja
				Chinese	一般强制操作标志
				Czech	Obecná značka pro povinnou akci
				Danish	Generelt skilt om obligatoriske handlinger
				Dutch	Algemeen verplichte actie-symbool
				Estonian	Üldine kohustusliku tegevuse märk
				Finnish	Yleinen pakollisen toiminnan merkki
				French	Symbole général de mesure obligatoire
				German	Allgemeines Gebotszeichen
				Greek	Γενική σήμανση υποχρεωτικής ενέργειας
				Hungarian	Általános kötelező intézkedés jele
				Italian	Segnale di azione generale obbligatoria
				Japanese	一般的な必須行動表示
				Latvian	Vispārīga obligātā darbība
				Lithuanian	Bendrasis privalomojo veiksmo ženklas
				Norwegian	Skilt for generelle obligatoriske tiltak
				Polish	Ogólny znak działania obowiązkowego
				Portuguese	Sinal geral de ação obrigatória
Romanian	Semn general de acțiune obligatorie				
Slovak	Symbol všeobecného povinného opatrenia				
Slovenian	Splošni znak za obveznost				
Spanish	Señal de obligación general				
Swedish	Allmän symbol för obligatorisk åtgärd				
Turkish	Genel zorunlu eylem işareti				


IEC 60417


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 <p>Example</p>  <p>5V 0.1A</p>	Two combined symbols 5034 (input), and 5031 (direct current)	ISO 7000 / IEC Graphical Symbols for use on equipment	Direct current input	English	Direct current input
				Bulgarian	Вход за постоянен ток
				Croatian	Ulaz istosmjernje struje
				Chinese	直流输入
				Czech	Vstup stejnosměrného proudu
				Danish	Jævnstrømsindgang
				Dutch	Gelijkstroomingang
				Estonian	Alalisvoolu sisend
				Finnish	Tasavirran tulo
				French	Entrée de courant continu
				German	Gleichstromeingang
				Greek	Είσοδος συνεχούς ρεύματος
				Hungarian	Egyenáramú bemenet
				Italian	Ingresso corrente continua
				Japanese	直流入力
				Latvian	Līdzstrāvas ievade
				Lithuanian	Nuolatinės srovės įvestis
				Norwegian	Likestrømsinngang
				Polish	Wejście prądu stałego
				Portuguese	Entrada de corrente contínua
Romanian	Intrare curent continuu				
Slovak	Vstup jednosmerného prúdu				
Slovenian	Enosmerni vhodni tok				
Spanish	Entrada de corriente continua				
Swedish	Likström in				
Turkish	Düz akım girişi				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	Two combined symbols 5035 (output), and 5031 (direct current)	ISO 7000 / IEC Graphical Symbols for use on equipment	Direct current output	English	Direct current output
				Bulgarian	Изход за постоянен ток
				Croatian	Izlaz istosmjerne struje
				Chinese	直流输出
				Czech	Výstup stejnosměrného proudu
				Danish	Jævnstrømsudgang
				Dutch	Gelijkstroomuitgang
				Estonian	Alalisvoolu väljund
				Finnish	Tasavirran lähtö
				French	Sortie de courant continu
				German	Gleichstromausgang
				Greek	Έξοδος συνεχούς ρεύματος
				Hungarian	Egyenáramú kimenet
				Italian	Uscita corrente continua
				Japanese	直流出力
				Latvian	Līdzstrāvas izvade
				Lithuanian	Nuolatinės srovės išvestis
				Norwegian	Likestrømsutgang
				Polish	Wyjście prądu stałego
				Portuguese	Saída de corrente contínua
Romanian	Ieșire curent continuu				
Slovak	Výstup jednosmerného prúdu				
Slovenian	Enosmerni izhodni tok				
Spanish	Salida de corriente continua				
Swedish	Likström ut				
Turkish	Düz akım çıkışı				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5019	ISO 7000 / IEC Graphical Symbols for use on equipment	Identifies a protective ground connector	English	Protective earth (ground)
				Bulgarian	Защитно заземяване (земя)
				Croatian	Zaštitni vodič (uzemljenje)
				Chinese	保护接地 (接地)
				Czech	Ochranné uzemnění (kostra)
				Danish	Beskyttende jordforbindelse
				Dutch	Veiligheidsaarding
				Estonian	Kaitsemaandus
				Finnish	Suojamaadoitus
				French	Mise à la terre protectrice (masse)
				German	Schutzleiter (Masse)
				Greek	Προστατευτική γείωση
				Hungarian	Védőföldelés
				Italian	Massa di protezione (terra)
				Japanese	保護接地 (アース)
				Latvian	Aizsargzemējums
				Lithuanian	Apsauginis įžeminimas (žemė)
				Norwegian	Jordingstilkobling
				Polish	Uziemienie ochronne
				Portuguese	Ligação terra de proteção
				Romanian	Împământare de protecție (legare la masă)
Slovak	Ochranné uzemnenie (zem)				
Slovenian	Ozemljitev (zemlja)				
Spanish	Puesta a tierra de protección (masa)				
Swedish	Skyddsjord				
Turkish	Koruyucu topraklama (toprak)				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5021	ISO 7000 / IEC Graphical Symbols for use on equipment	Identifies the terminals that bring the various parts of a system to the same potential	English	Equipotentiality
				Bulgarian	ЕКВИПОТЕНЦИАЛНОСТ
				Croatian	Izjednačavanje potencijala
				Chinese	等位性
				Czech	Ekvipotenciální uzemnění
				Danish	Ækvipotentialitet
				Dutch	Potentiaalvereffening
				Estonian	Potentsiaalühtlustus
				Finnish	Potentiaalintasaus
				French	Équipotentialité
				German	Potenzialausgleich
				Greek	Ισοδυναμικότητα
				Hungarian	Potenciálkiegyenlítés
				Italian	Equipotenzialità
				Japanese	等電位性
				Latvian	Sprieguma izlīdzināšana
				Lithuanian	Ekvipotencialumas
				Norwegian	Potensialutligning
				Polish	Wyrównanie potencjałów
				Portuguese	Equipotencialidade
				Romanian	Echipotențializare
Slovak	Vyrovnanie potenciálov				
Slovenian	Ekvipotencialnost				
Spanish	Equipotencialidad				
Swedish	Potentialutjämning				
Turkish	Eşpotansiyellik				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5032	ISO 7000 / IEC Graphical Symbols for use on equipment	Indicates that the equipment is for alternating current only	English	Alternating current
				Bulgarian	Променлив ток
				Croatian	Izmjenična struja
				Chinese	交流
				Czech	Střídavý proud
				Danish	Vekselstrøm
				Dutch	Wisselstroom
				Estonian	Vahelduvvool
				Finnish	Vaihtovirta
				French	Courant alternatif
				German	Wechselstrom
				Greek	Εναλλασσόμενο ρεύμα
				Hungarian	Váltóáram
				Italian	Corrente alternata
				Japanese	交流
				Latvian	Mainstrāva
				Lithuanian	Kintamoji srovė
				Norwegian	Vekselstrøm
				Polish	Prąd przemienny
				Portuguese	Corrente alternada
Romanian	Curent alternativ				
Slovak	Striedavý prúd				
Slovenian	Izmenični tok				
Spanish	Corriente alterna				
Swedish	Växelström				
Turkish	Alternatif akım				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5031	ISO 7000 / IEC Graphical Symbols for use on equipment	Indicates that the equipment is for direct current only	English	Direct current
				Bulgarian	Постоянен ток
				Croatian	Istosmjerna struja
				Chinese	直流
				Czech	Stejnoseměrný proud
				Danish	Jævnstrøm
				Dutch	Gelijkstroom
				Estonian	Alalisvool
				Finnish	Tasavirta
				French	Courant continu
				German	Gleichstrom
				Greek	Συνεχές ρεύμα
				Hungarian	Egyenáram
				Italian	Corrente continua
				Japanese	直流
				Latvian	Līdzstrāva
				Lithuanian	Nuolatinė srovė
				Norwegian	Likestrøm
				Polish	Prąd stały
				Portuguese	Corrente contínua
Romanian	Curent continuu				
Slovak	Jednosmerný prúd				
Slovenian	Enosmerni tok				
Spanish	Corriente continua				
Swedish	Likström				
Turkish	Düz akım				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5957	ISO 7000 / IEC Graphical Symbols for use on equipment	Indicates that the <i>medical device</i> is only to be used indoors	English	Only for indoor use
				Bulgarian	Само за използване на закрито
				Croatian	Samo za upotrebu u zatvorenim prostorima
				Chinese	仅在室内使用
				Czech	Pouze pro použití v interiéru
				Danish	Kun til indendørs brug
				Dutch	Alleen voor gebruik binnenshuis
				Estonian	Ainult siseruumides kasutamiseks
				Finnish	Vain sisäkäyttöön
				French	Réservé à un usage intérieur
				German	Nur zur Verwendung in Innenräumen
				Greek	Μόνο για χρήση σε εσωτερικούς χώρους
				Hungarian	Csak beltéri használatra
				Italian	Solo per uso in ambienti interni
				Japanese	室内使用向け
				Latvian	Lietošanai vienīgi telpās
				Lithuanian	Naudoti tik patalpoje
				Norwegian	Kun for innendørs bruk
				Polish	Wyłącznie do użytku w budynkach
				Portuguese	Apenas para utilização no interior
Romanian	Exclusiv pentru utilizare în interior				
Slovak	Len na použitie v interiéri				
Slovenian	Samo za notranjo uporabo				
Spanish	Solo para uso en interiores				
Swedish	Endast avsedd för inomhusbruk				
Turkish	Sadece kapalı alanlarda kullanım içindir				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5172	ISO 7000 / IEC Graphical Symbols for use on equipment	Indicates that the <i>medical device</i> is a Class II equipment	English	Class II equipment
				Bulgarian	Оборудване от клас II
				Croatian	Oprema klase II
				Chinese	二类设备
				Czech	Zařízení třídy II
				Danish	Klasse II-udstyr
				Dutch	Klasse II-apparatuur
				Estonian	II klassi seade
				Finnish	Luokan II laite
				French	Équipement de classe II
				German	Gerät der Klasse II
				Greek	Εξοπλισμός Τάξης II
				Hungarian	II. osztályba tartozó berendezés
				Italian	Apparecchiatura di classe II
				Japanese	クラス II 機器
				Latvian	II klases aprīkojums
				Lithuanian	II klasės įranga
				Norwegian	Klasse II-utstyr
				Polish	Wyposażenie klasy II
				Portuguese	Equipamento da Classe II
				Romanian	Echipament de clasa II
				Slovak	Zariadenie triedy II
				Slovenian	Oprema razreda II
Spanish	Equipo de clase II				
Swedish	Klass II-utrustning				
Turkish	II. sınıf ekipman				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	5009	IEC 60417	Identifies the switch to be used to turn the <i>medical device</i> in the stand-by position	English	Stand-by
				Bulgarian	РЕЖИМ НА ГОТОВНОСТ
				Croatian	Stanje pripravnosti
				Chinese	待机
				Czech	Pohotovostní režim
				Danish	Standby
				Dutch	Stand-by
				Estonian	Ooterežiim
				Finnish	Valmiustila
				French	En attente
				German	Stand-by
				Greek	Αναμονή
				Hungarian	Készlet
				Italian	Standby
				Japanese	スタンバイ
				Latvian	Gaidstāve
				Lithuanian	Budėjimo režimas
				Norwegian	Ventemodus
				Polish	Gotowość
				Portuguese	Em espera
Romanian	Standby				
Slovak	Pohotovostný režim				
Slovenian	Stanje pripravljenosti				
Spanish	Modo de espera				
Swedish	Standby				
Turkish	Beklemede				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 <ul style="list-style-type: none"> Green indicates that the power is On Orange indicates that the device is on stand-by mode 	5009	IEC/EN 60601-1	<p>A control device or switch that brings the ME EQUIPMENT into the “stand-by” condition may be indicated by use of symbol IEC 60417-5009</p> <p>The colors of indicator lights and their meanings shall comply with</p> <ul style="list-style-type: none"> Red: Warning immediate response by the OPERATOR is required Yellow: Caution prompt response by the OPERATOR is required Green: Ready for use <p>Any other color meaning other than that of red, yellow, green, or cyan.</p>	English	Status indicator
				Bulgarian	Индикатор за състоянието
				Croatian	Pokazivač stanja
				Chinese	状态指示灯
				Czech	Indikátor stavu
				Danish	Statusindikator
				Dutch	Statusindicator
				Estonian	Olekunäidik
				Finnish	Tilan ilmaisin
				French	Indicateur d'état
				German	Statusanzeige
				Greek	Δείκτης κατάστασης
				Hungarian	Állapotjelző
				Italian	Indicatore di stato
				Japanese	ステータスインジケータ
				Latvian	Statusa indikators
				Lithuanian	Būsenos indikatorius
				Norwegian	Status-indikatorer
				Polish	Wskaźnik stanu
				Portuguese	Indicador de estado
Romanian	Indicator de stare				
Slovak	Indikátor stavu				
Slovenian	Indikator stanja				
Spanish	Indicador de estado				
Swedish	Statusindikator				
Turkish	Durum göstergesi				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	4.3.8	IEC 60950:2013	N/A	English	Polarity of DC power connector
				Bulgarian	Полярност на конектор за постояннотоково захранване
				Croatian	Polaritet priključka za napajanje istosmjernom strujom
				Chinese	直流电源接头的极性
				Czech	Polarita konektoru stejnosměrného napájení
				Danish	Jævnstrømskonnektorens polaritet
				Dutch	Polariteit van DC-voedingsconnector
				Estonian	Alalisvoolupistikute polaarsus
				Finnish	DC-virtaliittimen napaisuus
				French	Polarité du connecteur d'alimentation CC
				German	Polarität des Gleichstromanschlusses
				Greek	Πολικότητα συνδέσμου ισχύος DC
				Hungarian	A DC-tápcsatlakozó polaritása
				Italian	Polarità del connettore di alimentazione CC
				Japanese	直流電源コネクタの極性
				Latvian	Līdzstrāvas savienotāja polaritāte
				Lithuanian	Nuolatinės srovės maitinimo jungties poliškumas
				Norwegian	Likestrømskontaktens polaritet
				Polish	Biegunowość złącza zasilania prądem stałym
				Portuguese	Polaridade do conector de alimentação CC
Romanian	Polaritatea conectorului de alimentare cu c.c.				
Slovak	Polarita DC napájacieho konektora				
Slovenian	Polarnost priključka enosmerne napetosti				
Spanish	Polaridad del conector de alimentación de CC				
Swedish	Likströmsanslutningens polaritet				
Turkish	DC güç konektörünün polaritesi				


UN3481 and UN3480


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	UN legislation - UN3481, Lithium-ion batteries packed with or contained in equipment	Indicates that the equipment is packed with or contains lithium-ion batteries	English	Lithium-ion batteries packed with or contained in equipment
				Bulgarian	Литиевойонни батерии, приложени към оборудването или съдържащи се в него
				Croatian	Litij-ionske baterije pakirane uz opremu ili sadržane u njoj
				Chinese	锂离子电池与设备包装在一起或内含于设备中
				Czech	Lithium-iontové baterie jako součást balení nebo obsažené v zařízení
				Danish	Litium-ion-batterier pakket med eller indeholdt i udstyret
				Dutch	Lithium-ionbatterijen verpakt met of omvat in apparatuur
				Estonian	Seadmega on kaasas või seadmesse on paigaldatud liitium-ioonakud
				Finnish	Pakkauksessa olevat tai laitteen sisältämät litiumioniakut
				French	Batteries lithium-ion emballées avec ou contenues dans l'équipement
				German	Im Gerät enthaltene oder mitgelieferte Lithium-Ionen-Akkus
				Greek	Η συσκευασία ή ο ίδιος ο εξοπλισμός περιέχουν μπαταρίες ιόντων λιθίου
				Hungarian	A berendezéshez mellékelte vagy benne található lítiumion-akkumulátorok
				Italian	Batterie agli ioni di litio confezionate o contenute nell'apparecchiatura
				Japanese	機器に内蔵または同梱されたリチウムイオン電池
				Latvian	Iekārtā vai tās iepakojumā ir litija jonu baterijas
				Lithuanian	Prie įrangos pridedamos arba joje yra ličio jonų baterijos
				Norwegian	Litium-ionbatterier pakket sammen med eller i utstyret
				Polish	Baterie litowo-jonowe do urządzenia lub zawarte w nim
				Portuguese	Baterias de íões de lítio embaladas com ou incluídas no equipamento
Romanian	Baterii litiu-ion ambalate cu sau conținute în echipament				
Slovak	Lítium-iónové batérie dodávané spolu so zariadením alebo obsiahnuté v zariadení				
Slovenian	Litij-ionske baterije, pakirane z opremo ali vsebovane v njej				


				Spanish	Baterías de iones de litio suministradas con el equipo o incluidas en él
				Swedish	Litiumjonbatterier förpackade med eller placerade i utrustningen
				Turkish	Ekipmanla birlikte paketlenmiş veya ekipmanın içinde bulunan lityum-iyon piller


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	UN legislation - UN3480, Lithium-ion batteries	Indicates lithium-ion batteries	English	Lithium-ion batteries
				Bulgarian	Литиевойонни батерии
				Croatian	Litij-ionske baterije
				Chinese	锂离子电池
				Czech	Lithium-iontové baterie
				Danish	Litium-ion-batterier
				Dutch	Lithium-ionbatterijen
				Estonian	Liitiumioonakud
				Finnish	Litiumioniakut
				French	Batteries lithium-ion
				German	Lithium-Ionen-Akkus
				Greek	Μπαταρίες ιόντων λιθίου
				Hungarian	Lítiumion-akkumulátorok
				Italian	Batterie agli ioni di litio
				Japanese	リチウムイオン電池
				Latvian	Litija jonu baterijas
				Lithuanian	Ličio jonų baterijos
				Norwegian	Litium-ionbatterier
				Polish	Baterie litowo-jonowe
				Portuguese	Baterias de iões de lítio
Romanian	Baterie litiu-ion				
Slovak	Lítium-iónové batérie				
Slovenian	Litij-ionske baterije				
Spanish	Baterías de iones de litio				
Swedish	Litiumjonbatterier				
Turkish	Lityum-iyon piller				


Symbols designed by AMBU


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that moving parts are not to be touched	English	Do not touch moving parts
				Bulgarian	Не пипайте движещи се части
				Croatian	Nemojte dodirivati pokretne dijelove
				Chinese	请勿触摸移动部件
				Czech	Nedotýkejte se pohyblivých částí
				Danish	Berør ikke bevægelige dele
				Dutch	Raak bewegende delen niet aan
				Estonian	Liikuvaid osi mitte puudutada
				Finnish	Älä koske liikkuviin osiin
				French	Ne pas toucher les pièces mobiles
				German	Bewegliche Teile nicht berühren
				Greek	Μην αγγίζετε τα κινούμενα μέρη
				Hungarian	Ne érintse meg a mozgó alkatrészeket!
				Italian	Non toccare le parti mobili
				Japanese	可動部品に触れないでください
				Latvian	Nepieskarties kustīgajām daļām
				Lithuanian	Nelieskite judančių dalių
				Norwegian	Ikke berør bevegelige deler
				Polish	Nie dotykać ruchomych części
				Portuguese	Não toque em peças em movimento
Romanian	Nu atingeți componentele în mișcare.				
Slovak	Nedotýkajte sa pohyblivých súčastí				
Slovenian	Ne dotikajte se premikajočih delov				
Spanish	No toque las piezas móviles				
Swedish	Vidrör inte rörliga delar				
Turkish	Hareketli parçalara dokunmayın				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the item must not be used more than 30 times	English	Do not use more than 30 times
				Bulgarian	Не използвайте повече от 30 пъти
				Croatian	Nemojte upotrebljavati više od 30 puta
				Chinese	使用次数不得超过 30 次
				Czech	Nepoužívejte více než 30krát
				Danish	Må ikke anvendes mere end 30 gange
				Dutch	Niet vaker dan 30 keer gebruiken
				Estonian	Mitte kasutada üle 30 korra
				Finnish	Käyttö sallittu enintään 30 kertaa
				French	Ne pas utiliser plus de 30 fois
				German	Nicht mehr als 30-mal verwenden
				Greek	Μην το χρησιμοποιείτε περισσότερες από 30 φορές
				Hungarian	Ne használja 30-nál több alkalommal!
				Italian	Non utilizzare più di 30 volte
				Japanese	30 回以上使用しないでください
				Latvian	Nelietot vairāk par 30 reizēm
				Lithuanian	Nenaudoti daugiau nei 30 kartų
				Norwegian	Ikke bruk mer enn 30 ganger
				Polish	Nie używać więcej niż 30 razy
				Portuguese	Não utilizar mais de 30 vezes
				Romanian	A nu se utiliza de mai mult de 30 ori
Slovak	Nepoužívajte viac ako 30-krát				
Slovenian	Ne uporabljajte več kot 30-krat				
Spanish	No utilizar más de 30 veces				
Swedish	Får användas högst 30 gånger				
Turkish	En fazla 30 kez kullanın				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the item must not be used more than 40 times	English	Do not use more than 40 times
				Bulgarian	Не използвайте повече от 40 пъти
				Croatian	Nemojte upotrebljavati više od 40 puta
				Chinese	使用次数不得超过 40 次
				Czech	Nepoužívejte více než 40krát
				Danish	Må ikke anvendes mere end 40 gange
				Dutch	Niet vaker dan 40 keer gebruiken
				Estonian	Mitte kasutada üle 40 korra
				Finnish	Käyttö sallittu enintään 40 kertaa
				French	Ne pas utiliser plus de 40 fois
				German	Nicht mehr als 40-mal verwenden
				Greek	Μην το χρησιμοποιείτε περισσότερες από 40 φορές
				Hungarian	Ne használja 40-nél több alkalommal!
				Italian	Non utilizzare più di 40 volte
				Japanese	40 回以上使用しないでください
				Latvian	Nelietot vairāk par 40 reizēm
				Lithuanian	Nenaudoti daugiau nei 40 kartų
				Norwegian	Ikke bruk mer enn 40 ganger
				Polish	Nie używać więcej niż 40 razy
				Portuguese	Não utilizar mais de 40 vezes
				Romanian	A nu se utiliza de mai mult de 40 ori
Slovak	Nepoužívajte viac ako 40-krát				
Slovenian	Ne uporabljajte več kot 40-krat				
Spanish	No utilizar más de 40 veces				
Swedish	Får användas högst 40 gånger				
Turkish	En fazla 40 kez kullanın				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the item must not be used more than 100 times	English	Do not use more than 100 times
				Bulgarian	Не използвайте повече от 100 пъти
				Croatian	Nemojte upotrebljavati više od 100 puta
				Chinese	使用次数不得超过 100 次
				Czech	Nepoužívejte více než 100krát
				Danish	Må ikke anvendes mere end 100 gange
				Dutch	Niet vaker dan 100 keer gebruiken
				Estonian	Mitte kasutada üle 100 korra
				Finnish	Käyttö sallittu enintään 100 kertaa
				French	Ne pas utiliser plus de 100 fois
				German	Nicht mehr als 100-mal verwenden
				Greek	Μην το χρησιμοποιείτε περισσότερες από 100 φορές
				Hungarian	Ne használja 100-nél több alkalommal!
				Italian	Non utilizzare più di 100 volte
				Japanese	100 回以上使用しないでください
				Latvian	Nelietot vairāk par 100 reizēm
				Lithuanian	Nenaudoti daugiau nei 100 kartų
				Norwegian	Ikke bruk mer enn 100 ganger
				Polish	Nie używać więcej niż 100 razy
				Portuguese	Não utilizar mais de 100 vezes
				Romanian	A nu se utiliza de mai mult de 100 ori
Slovak	Nepoužívajte viac ako 100-krát				
Slovenian	Ne uporabljajte več kot 100-krat				
Spanish	No utilizar más de 100 veces				
Swedish	Får användas högst 100 gånger				
Turkish	En fazla 100 kez kullanın				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the pouch must be kept closed using a clip	English	Open pouch must be closed with a clip
				Bulgarian	Отворената торбичка трябва да бъде затворена с клипс
				Croatian	Otvorena vrećica mora se zatvoriti korčom
				Chinese	袋子打开后必须用夹子封好
				Czech	Otevřený vak musí být uzavřen svorkou
				Danish	Den åbne pose skal lukkes med en clips
				Dutch	Open zakje moet worden gesloten met een clip
				Estonian	Avatud kott tuleb klambriga sulgeda
				Finnish	Avoin pussi on suljettava sulkimella
				French	Le sachet ouvert doit être fermé à l'aide d'une pince
				German	Offener Beutel muss mit einem Clip verschlossen werden
				Greek	Η ανοικτή θήκη πρέπει να κλείνει με κλιπ
				Hungarian	A nyitott tasakot kapoccsal le kell zárni
				Italian	Chiudere la busta aperta con una clip
				Japanese	オープンパウチは、クリップで閉じる必要があります
				Latvian	Atvērtais maisiņš ir jāaizver ar klipsi
				Lithuanian	Atidarytas maišelis turi būti uždarytas spaustuku
				Norwegian	Åpnet pose må lukkes med en klemme
				Polish	Otwarty woreczek musi być zamknięty klamrą
				Portuguese	A bolsa aberta tem de estar fechada com uma mola
				Romanian	Punga deschisă trebuie închisă cu o clemă
				Slovak	Otvorené vrečko sa musí uzavrieť svorkou
				Slovenian	Odprto vrečko je treba zapreti s sponko
Spanish	La bolsa abierta debe cerrarse con un clip				
Swedish	Öppen påse måste förslutas med en klämma				
Turkish	Açık torba, bir klipsle kapatılmalıdır				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the <i>medical device</i> permits the x-ray beam to pass through it	English	X-ray Radiolucent
				Bulgarian	Рентгенопрозрачно
				Croatian	Rendgenski radiolucentno
				Chinese	可透 X 射线
				Czech	Prostupné pro RTG záření
				Danish	Røntgentransparent
				Dutch	Röntgendoorlaatbaar
				Estonian	Röntgenläbipaistev
				Finnish	Röntgennegatiivinen
				French	Radiotransparent aux rayons X
				German	Röntgendurchlässig
				Greek	Ακτινοδιαπερατό στις ακτίνες X
				Hungarian	Áteresztí a röntgensugarakat
				Italian	Radiotrasparente ai raggi X
				Japanese	X 線透過
				Latvian	Rentgenstarus caurlaidīgs
				Lithuanian	Pralaidus rentgeno spinduliams
				Norwegian	Røntgentransparent
				Polish	Przepuszcza promienie rentgenowskie
				Portuguese	Raio-X radiolucente
				Romanian	Radiotransparent la raze X
Slovak	Prepúšťa röntgenové žiarenie				
Slovenian	Prepušča rentgenske žarke				
Spanish	Radiografía radiotransparente				
Swedish	Röntgengenomsläpplig				
Turkish	Radyolusent X ışını				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the <i>medical device</i> contains solid materials	English	Solid
				Bulgarian	Твърдо
				Croatian	Čvrsto
				Chinese	固体
				Czech	Pevný materiál
				Danish	Fastgel
				Dutch	Vast
				Estonian	Tahke
				Finnish	Kiinteä
				French	Solide
				German	Fest
				Greek	Στερεό
				Hungarian	Szilárd
				Italian	Solido
				Japanese	ソリッド
				Latvian	Ciets
				Lithuanian	Kietas
				Norwegian	Fast
				Polish	Ciała stałe
				Portuguese	Sólido
Romanian	Solid				
Slovak	Tuhé				
Slovenian	Trdno				
Spanish	Sólido				
Swedish	Fast				
Turkish	Katı				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the <i>medical device</i> contains cloth materials	English	Cloth backing
				Bulgarian	Текстилна основа
				Croatian	Podloga od tkanine
				Chinese	布面基底
				Czech	Textilní podklad
				Danish	Tekstilbagside
				Dutch	Achterlaag van stof
				Estonian	Riidest kihiga
				Finnish	Kangas
				French	Support en tissu
				German	Vlies-Trägermaterial
				Greek	Υφασμάτινο υπόστρωμα
				Hungarian	Textil anyagot tartalmaz
				Italian	Rivestimento in tessuto
				Japanese	クロス裏当て
				Latvian	Auduma oderējums
				Lithuanian	Audinio pamušalas
				Norwegian	Bakside av tekstil
				Polish	Podkład materiałowy
				Portuguese	Suporte em tecido
Romanian	Căptușeală din pânză				
Slovak	Tkaninový podklad				
Slovenian	Tekstilna obloga				
Spanish	Respaldo de tejido				
Swedish	Textilstoppning				
Turkish	Kumaş arka destek				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the <i>medical device</i> contains microporous materials	English	Microporous backing
				Bulgarian	Микропореста основа
				Croatian	Mikroporozna podloga
				Chinese	微孔基底
				Czech	Mikroporézni podklad
				Danish	Mikroporøs bagside
				Dutch	Microporeuze achterlaag
				Estonian	Mikropoorse kihiga
				Finnish	Kuitukangas
				French	Support microporeux
				German	Mikroporöses Trägermaterial
				Greek	Υπόστρωμα από μικροπορώδες υλικό
				Hungarian	Mikropórusos anyagot tartalmaz
				Italian	Rivestimento microporoso
				Japanese	微孔性裏当て
				Latvian	Mikroporains oderējums
				Lithuanian	Mikroporinės medžiagos pamušalas
				Norwegian	Mikroporøs bakside
				Polish	Podkład mikroporowaty
				Portuguese	Suporte microporoso
				Romanian	Căptușeală din material microporos
Slovak	Mikroporézny podklad				
Slovenian	Mikroporna obloga				
Spanish	Respaldo microporoso				
Swedish	Mikroporøs stoppning				
Turkish	Mikro gözenekli arka destek				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the <i>medical device</i> contains wet materials	English	Wet
				Bulgarian	Мокро
				Croatian	Mokro
				Chinese	湿
				Czech	Vlhký materiál
				Danish	Våd
				Dutch	Nat
				Estonian	Märg
				Finnish	Märkä
				French	Liquide
				German	Nassgel
				Greek	Υγρό
				Hungarian	Nedves
				Italian	Liquido
				Japanese	ウェット
				Latvian	Mitrs
				Lithuanian	Drėgnas
				Norwegian	Våt
				Polish	Mokre
				Portuguese	Húmido
				Romanian	Ud
Slovak	Mokrý				
Slovenian	Mokro				
Spanish	Líquido				
Swedish	Våt				
Turkish	Islak				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the <i>medical device</i> contains foam materials	English	Foam backing
				Bulgarian	Пенопластова основа
				Croatian	Podloga od pjene
				Chinese	泡沫基底
				Czech	Pěnový podklad
				Danish	Skumbagside
				Dutch	Achterlaag van schuim
				Estonian	Vahtmaterjalist kihiga
				Finnish	Vaahтомуovi
				French	Support en mousse
				German	Schaumstoffträger
				Greek	Υπόστρωμα από αφρώδες υλικό
				Hungarian	Hab anyagot tartalmaz
				Italian	Supporto in foam
				Japanese	フォーム裏当て
				Latvian	Putuplasta oderējums
				Lithuanian	Putplasčio pamušalas
				Norwegian	Bakside av skum
				Polish	Podkład piankowy
				Portuguese	Suporte em espuma
				Romanian	Căptușeală din spumă
Slovak	Penový podklad				
Slovenian	Penasta obloga				
Spanish	Respaldo de espuma				
Swedish	Skumstopning				
Turkish	Köpük arka destek				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 Max OD	N/A	N/A	Indicates the outer diameter of the <i>medical device's</i> insertion portion	English	Maximum insertion portion width (Maximum outer diameter)
				Bulgarian	Максимална ширина на въвежданата част (максимален външен диаметър)
				Croatian	Maksimalna širina uvodnog dijela (maksimalni vanjski promjer)
				Chinese	插入部最大宽度 (最大外径)
				Czech	Maximální šířka zavedené části (maximální vnější průměr)
				Danish	Maksimal bredde på indført del (maksimal ydre diameter)
				Dutch	Maximale breedte inbrengdeel (maximale buitendiameter)
				Estonian	Sisestatava osa maksimaalne laius (maksimaalne väline diameeter)
				Finnish	Sisäänviedyn osan enimmäisleveys (maksimiulkohalkaisija)
				French	Largeur maximale de la partie insérée (diamètre extérieur maximal)
				German	Maximale Breite des Intubationsbereichs (maximaler Außendurchmesser)
				Greek	Μέγιστο πλάτος τμήματος εισαγωγής (μέγιστη εξωτερική διάμετρος)
				Hungarian	A behelyezett rész maximális szélessége (max. külső átmérő)
				Italian	Larghezza massima della parte d'inserzione (diametro esterno massimo)
				Japanese	最大挿入部幅 (最大外径)
				Latvian	Maksimālais ievadišanas daļas platumš (maksimālais ārējais diametrs)
				Lithuanian	Didžiausias įvedamos dalies plotis (maksimalus išorinis skersmuo)
				Norwegian	Maksimal bredde på innført del (maksimal utvendig diameter)
				Polish	Maksymalna szerokość wprowadzanego odcinka (maksymalna średnica zewnętrzna)
				Portuguese	Largura máxima da parte inserida (diâmetro externo máximo)
				Romanian	Lățime maximă a porțiunii de introduce (diametru exterior maxim)
				Slovak	Maximálna šírka zavádzacej časti (maximálny vonkajší priemer)
				Slovenian	Največja širina vstavljenega dela (največji zunanji premer)
Spanish	Ancho máximo de la parte insertada (diámetro exterior máximo)				
Swedish	Maximal bredd för införingsdel (max. ytterdiameter)				
Turkish	Yerleşirme kısmının maksimum genişliği (Maksimum dış çap)				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 Min ID	N/A	N/A	Indicates the inner diameter of the <i>medical device's</i> insertion portion	English	Minimum instrument channel width (Minimum inner diameter)
				Bulgarian	Минимална ширина на канала на инструмента (минимален вътрешен диаметър)
				Croatian	Minimalna širina kanala instrumenta (minimalni unutarnji promjer)
				Chinese	仪器通道最小宽度 (最小内径)
				Czech	Minimální šířka nástrojového kanálu (minimální vnitřní průměr)
				Danish	Minimumsbredde på instrumentkanal (minimal indvendig diameter)
				Dutch	Minimale breedte van het instrumentkanaal (minimale binnendiameter)
				Estonian	Instrumendi kanali minimaalne laius (minimaalne sisemine diameeter)
				Finnish	Instrumentikanavan vähimmäisleveys (minimisisähalkaisija)
				French	Largeur minimale du canal pour instruments (diamètre intérieur minimal)
				German	Minimale Breite des Instrumentenkanals (minimaler Innendurchmesser)
				Greek	Ελάχιστο πλάτος καναλιού οργάνου (ελάχιστη εσωτερική διάμετρος)
				Hungarian	A műszercsatorna minimális szélessége (min. belső átmérő)
				Italian	Larghezza minima del canale dello strumento (diametro interno minimo)
				Japanese	最小機器チャンネル幅 (最小内径)
				Latvian	Instrumentu kanāla minimālais platums (minimālais iekšējais diametrs)
				Lithuanian	Mažiausias instrumento kanalo plotis (minimalus vidinis skersmuo)
				Norwegian	Minimum bredde på instrumentkanal (minimum innvendig diameter)
				Polish	Minimalna szerokość kanału instrumentalnego (minimalna średnica wewnętrzna)
				Portuguese	Largura mínima do canal do instrumento (diâmetro interno mínimo)
Romanian	Lățimea minimă a canalului instrumentului (diametru interior minim)				
Slovak	Minimálna šírka kanála nástroja (minimálny vnútorný priemer)				
Slovenian	Najmanjša širina delovnega instrumenta (najmanjši notranji premer)				
Spanish	Ancho mínimo del canal del instrumento (diámetro interior mínimo)				
Swedish	Minimibredd för instrumentkanalen (minsta innerdiameter)				
Turkish	Alet kanalının minimum genişliği (Minimum iç çap)				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	FOV = 90° (Field of View)	English	Field of view
			FOV = 120° (Field of View)	Bulgarian	Зрително поле
				Croatian	Polje prikaza
				Chinese	视场
				Czech	Zorné pole
				Danish	Synsfelt
				Dutch	Gezichtsveld
				Estonian	Vaateväli
				Finnish	Kuvakulma
				French	Champ de vision
				German	Sichtfeld
				Greek	Οπτικό πεδίο
				Hungarian	Látómező
				Italian	Campo visivo
				Japanese	視野角
				Latvian	Skata lauks
				Lithuanian	Matomumo laukas
				Norwegian	Synsfelt
				Polish	Pole widzenia
				Portuguese	Campo de visão
				Romanian	Câmp de vizualizare
				Slovak	Zorné pole
				Slovenian	Zorno polje
				Spanish	Campo de visión
				Swedish	Betraktningsfält
				Turkish	Görüş alanı

Country-specific symbols
US and Canada

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the medical device has been evaluated by UL for a specific hazard or for performance under specific conditions	English	UL Classified Mark
				Bulgarian	Маркировка за класификация по UL
				Croatian	UL oznaka
				Chinese	UL 认证标志
				Czech	Značka klasifikace UL
				Danish	UL-mærket
				Dutch	Merkteken Door UL geclassificeerd
				Estonian	UL klassifikatsiooni märgis
				Finnish	UL-luokitusmerkintä
				French	Classification UL
				German	UL-Kennzeichnung
				Greek	Σήμανση ταξινόμησης κατά UL
				Hungarian	UL osztályozási jelölés
				Italian	Marchio classificazione UL
				Japanese	UL 分類マーク
				Latvian	UL klasifikācijas marķējums
				Lithuanian	UL klasifikacijos ženklas
				Norwegian	UL-klassifisert merke
				Polish	Znak klasyfikacji UL
				Portuguese	Marca de classificação UL
				Romanian	Marcaj UL
				Slovak	Značka klasifikácie podľa UL
Slovenian	Zaščitna oznaka UL				
Spanish	Marca de clasificación UL				
Swedish	UL-märkning				
Turkish	UL Sınıflandırılmış İşareti				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the component is part of a larger UL-certified product or system	English	UL Recognized Component Mark
				Bulgarian	Маркировка за признат компонент от UL
				Croatian	Oznaka UL priznate komponente
				Chinese	UL 认可组件标志
				Czech	Značka uznané součástky UL
				Danish	UL-anerkendt komponentmærke
				Dutch	Merkteken Door UL erkend onderdeel
				Estonian	UL tunnustatud komponendi märgis
				Finnish	UL-tunnustetun komponentin merkintä
				French	Marquage UL
				German	UL-Kennzeichnung für anerkannte Komponenten
				Greek	Σήμανση αναγνωρισμένων εξαρτημάτων κατά UL
				Hungarian	UL-minősítésű komponens jelölése
				Italian	Marchio componenti riconosciuti UL
				Japanese	UL 認定コンポーネントマーク
				Latvian	UL atzīta komponenta marķējums
				Lithuanian	UL pripažintas komponento ženklas
				Norwegian	UL-anerkjent komponentmerke
				Polish	Znak komponentu uznanego przez UL
				Portuguese	Marca de componente reconhecido pela UL
				Romanian	Marcaj UL pentru componente recunoscute
Slovak	Značka komponentu uznávaného UL				
Slovenian	Oznaka priznane komponente UL				
Spanish	Marca de componente reconocido por UL				
Swedish	Märkning för UL-godkända komponenter				
Turkish	UL Tanınan Bileşen İşareti				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the medical device has been tested and the electromagnetic radiation from the device is below the limits specified by the Federal Communications Commission (FCC).	English	Federal Communications Commission Mark
				Bulgarian	Маркировка на Федерална комисия по комуникациите
				Croatian	Oznaka Savezne komisije za komunikacije
				Chinese	联邦通信委员会标志
				Czech	Značka Federal Communications Commission
				Danish	Føderalt kommunikationskommissionsmærkning
				Dutch	Merkteken van de FCC
				Estonian	Föderaalse kommunikatsioonikomisjoni märgis
				Finnish	FCC-merkintä
				French	Marquage FCC
				German	Kennzeichen der Federal Communications Commission
				Greek	Σήμανση Ομοσπονδιακής Επιτροπής Επικοινωνιών
				Hungarian	Az USA Szövetségi Kommunikációs Bizottságának (FCC) jelölése
				Italian	Marchio Federal Communications Commission
				Japanese	連邦通信委員会マーク
				Latvian	Federālās sakaru komisijas marķējums
				Lithuanian	Federalinės ryšių komisijos ženklas
				Norwegian	Merke for Federal Communications Commission
				Polish	Znak Federalnej Komisji Łączności
				Portuguese	Marca da Comissão Federal de Comunicações
Romanian	Marcajul Comisiei Federale de Comunicații				
Slovak	Značka Federálnej komunikačnej komisie				
Slovenian	Oznaka Zvezne komisije za komunikacije				
Spanish	Marca de la Comisión Federal de Comunicaciones				
Swedish	Märke för Federal Communications Commission (FCC)				
Turkish	Federal İletişim Komisyonu İşareti				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
CONTAINS IC ID : 4104A-QCNFA324	N/A	N/A	The module is certified by Bureau Veritas, in relation to Canadian Radio Frequency requirements.	English	CAN: IC
				Bulgarian	CAN: IC
				Croatian	CAN: IC
				Chinese	CAN: IC
				Czech	KANADA: IC
				Danish	CAN: IC
				Dutch	CAN: IC
				Estonian	CAN : IC
				Finnish	CAN: IC
				French	CANADA : cause des interférences
				German	CAN: IC
				Greek	CAN: IC
				Hungarian	CAN: IC
				Italian	CAN: CI
				Japanese	CAN:IC
				Latvian	CAN: IC
				Lithuanian	CAN: IC
				Norwegian	CAN: IC
				Polish	CAN: IC
				Portuguese	CAN: IC
Romanian	CANADA: IC				
Slovak	KAN: Spôsobuje interferenciu				
Slovenian	CAN: IC				
Spanish	CAN: CI				
Swedish	CAN: IC				
Turkish	CAN: IC				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
CAN ICES-1/NMB-1	N/A	N/A	The radio frequency module has been Certified by: The module is certified by Bureau Veritas, in relation to Canadian Radio Frequency requirements.	English	N/A
				Bulgarian	Не е приложимо
				Croatian	Nije primjenjivo
				Chinese	不适用
				Czech	N/A
				Danish	Ikke tilgængelig
				Dutch	n.v.t.
				Estonian	Ei ole rakendatav
				Finnish	N/A
				French	N/A
				German	N/A
				Greek	Δ/Ι
				Hungarian	N/A
				Italian	N/A
				Japanese	該当なし
				Latvian	N/A
				Lithuanian	Netaikoma
				Norwegian	I/A
				Polish	Nd.
				Portuguese	N/A
Romanian	Nu este cazul				
Slovak	Neuvádza sa				
Slovenian	/				
Spanish	N/A				
Swedish	N/A				
Turkish	Yok				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
US: Rx only	(c) (1) (i) (F) (b) (1)	21 CFR 801.15 21 CFR 01.109	Caution: Federal (US) law restricts this device to sale by or on the order of a physician.	English	Caution: Federal law restricts this device to sale by or on the order of a licensed health care practitioner
				Bulgarian	Внимание: Федералното законодателство налага ограничението това изделие да се продава от или по поръчка на лицензиран практикуващ здравен специалист
				Croatian	Opresz: Savezni zakon ograničava prodaju ovog uređaja na prodaju od strane ili po nalogu ovlaštenog zdravstvenog djelatnika
				Chinese	注意事项：联邦法律限定，本器械只能由有执业资格的医疗从业人员销售，或者按照其医嘱购买
				Czech	Upozornění: Federální zákon (USA) omezuje prodej tohoto prostředku pouze lékařům nebo na předpis lékaře
				Danish	Bemærk: I henhold til amerikansk lovgivning må dette udstyr kun sælges af en autoriseret læge eller på foranledning af en autoriseret læge
				Dutch	Voorzichtig: Federale wetgeving staat de verkoop van dit apparaat uitsluitend toe door of op voorschrift van een arts of andere professional met een medische vergunning
				Estonian	Hoiatus! Föderaalseedused lubavad seda seadet müüa ainult litsentseeritud tervishoiutöötjate poolt või nende korraldusel
				Finnish	Huomio: Liittovaltion lainsäädäntö rajoittaa tämän laitteen myynnin ja tilaamisen asianmukaisesti valtuutetuille terveydenhuollon ammattilaisille
				French	Attention : les lois fédérales exigent que la vente de ce dispositif ne se fasse que sur ordonnance ou par l'intermédiaire d'un professionnel de santé autorisé
				German	Achtung: Gemäß Bundesrecht ist dieses Gerät beschränkt auf den Verkauf durch oder auf Anordnung eines Arztes oder einer ordnungsgemäß lizenzierten medizinischen Fachkraft
Greek	Προσοχή: Η ομοσπονδιακή νομοθεσία περιορίζει την πώληση του παρόντος ιατροτεχνολογικού προϊόντος από ιατρό ή κατόπιν εντολής ιατρού				


				Hungarian	Vigyázat! A szövetségi törvények értelmében ez az eszköz kizárólag engedéllyel rendelkező egészségügyi szakember által vagy rendelvényére értékesíthető
				Italian	Attenzione: la legge federale degli Stati Uniti limita la vendita del presente dispositivo ai soli operatori sanitari abilitati
				Japanese	注意：連邦法では、この機器の販売を、有資格医療従事者によるものか、有資格医療従事者の指示によるものに限定しています
				Latvian	UZMANĪBU! Saskaņā ar federālo likumu šo ierīci drīkst pārdot tikai licencēts veselības aprūpes speciālists vai pēc licencēta speciālista pasūtījuma
				Lithuanian	Dėmesio: pagal federalinius įstatymus šį prietaisą gali parduoti tik licencijuotas sveikatos priežiūros specialistas arba pagal jo nurodymą
				Norwegian	Forsiktig: Føderal lov begrenser salget av enheten til å kun være tillatt av eller etter ordre fra autorisert helsepersonell
				Polish	Uwaga: Przepisy federalne zezwalają na sprzedaż tego urządzenia wyłącznie przez uprawnionego pracownika służby zdrowia lub na jego zlecenie
				Portuguese	Atenção: Em conformidade com a lei federal, este dispositivo só pode ser vendido por um profissional de saúde licenciado ou mediante prescrição do mesmo
				Romanian	Atenție! Potrivit legilor federale din SUA, vânzarea acestui dispozitiv se efectuează doar prin sau la comanda unui medic licențiat
				Slovak	Upozornenie: Federálne zákony (USA) povoľujú predaj tejto pomôcky len lekárom alebo na objednávku lekára
				Slovenian	Pozor: Zvezni zakon omejuje prodajo te naprave pooblaščenim zdravstvenim delavcem ali po njihovem naročilu
				Spanish	Precaución: La legislación federal solo permite la venta de este dispositivo a través de profesionales sanitarios autorizados o bajo prescripción facultativa
				Swedish	Försiktighet: Enligt federal lagstiftning i USA får denna produkt endast säljas av läkare eller på läkares ordination
				Turkish	Dikkat: Federal yasa uyarınca bu cihaz, sadece lisanslı bir hekim tarafından ya da onun talimatı üzerine satılabilir


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	NRTL symbol for TÜV SÜD	English	NRTL Symbol
				Bulgarian	NRTL символ
				Croatian	Simbol NRTL
				Chinese	NRTL 符号
				Czech	Symbol NRTL
				Danish	NRTL-symbol
				Dutch	NRTL-symbool
				Estonian	Riiklikult tunnustatud katselabori (NRTL) sümbol
				Finnish	NRTL-symboli
				French	Symbole NRTL
				German	NRTL-Symbol
				Greek	Σύμβολο NRTL
				Hungarian	NRTL-szimbólum
				Italian	Simbolo NRTL
				Japanese	NRTL シンボル
				Latvian	NRTL simbols
				Lithuanian	NRTL simbolis
				Norwegian	NRTL-symbol
				Polish	Symbol NRTL
				Portuguese	Símbolo NRTL
Romanian	Simbol NRTL				
Slovak	Symbol NRTL				
Slovenian	Simbol NRTL				
Spanish	Símbolo NRTL				
Swedish	NRTL-symbol				
Turkish	NRTL Sembolü				

Australia and New Zealand


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	AS/NZS 3820: Essential safety requirements for electrical equipment	Indicates that the <i>medical device</i> complies with essential safety requirements for electrical equipment in Australia and New Zealand	English	Regulatory Compliance Mark
				Bulgarian	Маркировка за нормативно съответствие
				Croatian	Oznaka regulatorne usklađenosti
				Chinese	合规标志
				Czech	Značka shody s předpisy
				Danish	Regeloverensstemmelsesmærkning
				Dutch	Merkteken voor regelgevende naleving
				Estonian	Regulatiivse vastavuse märgis
				Finnish	RMC-merkintä, Australian ja Uuden-Seelannin lainsäädännön mukainen
				French	Marquage RCM
				German	Prüfzeichen für die Einhaltung gesetzlicher Vorschriften
				Greek	Σήμανση κανονιστικής συμμόρφωσης
				Hungarian	Megfelelőségi szabályozási jelölés
				Italian	Marchio di conformità normativa
				Japanese	規制コンプライアンスマーク
				Latvian	Normatīvo prasību ievērošanas marķējums
				Lithuanian	Reguliavimo atitikties ženklas
				Norwegian	Merke for regeletterlevelse
				Polish	Znak zgodności z przepisami
				Portuguese	Marca de conformidade regulamentar
				Romanian	Marcaj de conformitate cu reglementările
				Slovak	Značka zhody s predpismi
				Slovenian	Oznaka regulativne skladnosti
Spanish	Marca de conformidad reglamentaria				
Swedish	Regulatory Compliance Mark				
Turkish	Mevzuata Uygunluk İşareti				

Japan



Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the medical device complies with the requirements stated for Specified Radio Equipment in the Japanese Radio Law	English	Technical Standards Conformity Mark
				Bulgarian	Маркировка за съответствие с технически стандарти
				Croatian	Oznaka sukladnosti s tehničkim standardima
				Chinese	技术标准适合认定标志
				Czech	Značka shody s technickými normami
				Danish	Overensstemmelsesmærkning for teknisk standard
				Dutch	Merkteken van conformiteit met technische normen
				Estonian	Tehnilistele standarditele vastavuse märgis
				Finnish	Japanin teknisten standardien vaatimustenmukaisuusmerkintä
				French	Marquage Technical Standards Conformity
				German	Konformitätskennzeichen für technische Normen
				Greek	Σήμανση συμμόρφωσης κατά τα τεχνικά πρότυπα
				Hungarian	Műszaki szabványoknak való megfelelés jelölése
				Italian	Marchio conformità alle norme tecniche
				Japanese	技術規格適合性マーク
				Latvian	Tehnisko standartu atbilstības marķējums
				Lithuanian	Techninių standartų atitikties ženklas
				Norwegian	Merke for samsvar med tekniske standarder
				Polish	Znak zgodności z normami technicznymi
				Portuguese	Marca de conformidade com as normas técnicas
Romanian	Marcaj de conformitate cu standardele tehnice				
Slovak	Značka zhody s technickými normami				
Slovenian	Oznaka skladnosti s tehničnimi standardi				
Spanish	Marca de conformidad con las normas técnicas				
Swedish	Märkning för överensstämmelse med tekniska standarder				
Turkish	Teknik Standartlara Uygunluk İşareti				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the medical device complies with the electrical safety and EMI requirements for Japan	English	Japanese PSE marking
				Bulgarian	Японска маркировка PSE
				Croatian	Japanska oznaka PSE
				Chinese	日本 PSE 标志
				Czech	Japonská značka PSE
				Danish	Japansk PSE-mærkning
				Dutch	Japans PSE-markering
				Estonian	Jaapani PSE-märgis
				Finnish	Japanin PSE-merkintä
				French	Marquage PSE japonais
				German	Japanische PSE-Kennzeichnung
				Greek	Σήμανση ιαπωνικού PSE
				Hungarian	Japán PSE-jelölés
				Italian	Marchio PSE giapponese
				Japanese	日本 PSE マーク
				Latvian	PSE marķējums Japānā
				Lithuanian	Japonijos PSE ženklimas
				Norwegian	Japansk PSE-merking
				Polish	Japońskie oznakowanie PSE
				Portuguese	Marcação PSE japonesa
Romanian	Marcaj PSE pentru Japonia				
Slovak	Japonské označenie PSE				
Slovenian	Japonska oznaka PSE				
Spanish	Marcado PSE japonés				
Swedish	PSE-märkning, Japan				
Turkish	Japon PSE işaretlemeşi				


Taiwan


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates that the medical device complies with the minimum legal requirements for radio technology products	English	NCC Certification for Taiwan
				Bulgarian	NCC сертификат за Тайван
				Croatian	Certifikat NCC za Tajvan
				Chinese	台灣 NCC 认证
				Czech	Certifikace NCC pro Tchaj-wan
				Danish	NCC-certificering for Taiwan
				Dutch	NCC-certificering voor Taiwan
				Estonian	NCC-sertifikaat Taiwani jaoks
				Finnish	Taiwanin NCC-sertifiointi
				French	Certification NCC pour Taïwan
				German	NCC-Zertifizierung für Taiwan
				Greek	Πιστοποίηση NCC για Ταϊβάν
				Hungarian	Tajvanra vonatkozó NCC-tanúsítvány
				Italian	Certificazione NCC per Taiwan
				Japanese	台灣 NCC 認証
				Latvian	NCC sertifikāts Taivānā
				Lithuanian	NCC sertifikavimas Taivanui
				Norwegian	NCC-sertifisering for Taiwan
				Polish	Certyfikat NCC dla Tajwanu
				Portuguese	Certificação NCC para Taiwan
Romanian	Certificare NCC pentru Taiwan				
Slovak	Certifikácia NCC pre Taiwan				
Slovenian	Certificiranje NCC za Tajvan				
Spanish	Certificación NCC para Taiwán				
Swedish	NCC-certifiering, Taiwan				
Turkish	Tayvan için NCC Belgelendirmesi				

Brazil

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
<p>Segurança</p>  <p>Segurança</p> 	N/A	INMETRO 3.1 edition of IEC 60601-1	Indicates that the <i>medical device</i> meets Brazilian electrical safety requirements and is certified by an INMETRO accredited organization	English	INMETRO Certificate Medical Electrical Equipment
				Bulgarian	Сертификат INMETRO за електромедицински апарати
				Croatian	INMETRO certifikat za medicinsku električnu opremu
				Chinese	医用电气设备 INMETRO 认证
				Czech	Certifikát INMETRO pro zdravotnické elektrické přístroje
				Danish	INMETRO-certifikat for medicinsk elektrisk udstyr
				Dutch	INMETRO-certificaat inzake medische elektrische apparatuur
				Estonian	INMETRO sertifitseeritud meditsiiniseadmetik
				Finnish	INMETRO-sertifikaatti, lääkinnälliset sähkölaitteet
				French	Certificat INMETRO – Équipement électrique médical
				German	INMETRO-Zertifizierung für medizinische elektrische Geräte
				Greek	INMETRO Πιστοποιητικό Ηλεκτρικού Ιατρικού εξοπλισμού
				Hungarian	Gyógyászati villamos készülékek INMETRO tanúsítványa
				Italian	Certificato INMETRO per apparecchiature elettromedicali
				Japanese	INMETRO 証明書 (医療用電気機器)
				Latvian	INMETRO sertifikāts medicīnas elektroiekārtām
				Lithuanian	INMETRO medicininės elektros įrangos sertifikatas
				Norwegian	INMETRO-certifikat medisinsk elektrisk utstyr
				Polish	Certyfikat INMETRO Medyczne urządzenia elektryczne
				Portuguese	Certificação INMETRO de equipamento elétrico médico
Romanian	Certificat INMETRO pentru echipamente electrice medicale				
Slovak	Zdravotnícke elektrické zariadenia certifikované podľa INMETRO				
Slovenian	Certifikat INMETRO za medicinsko električno opremo				
Spanish	Certificado INMETRO para equipos electromédicos				
Swedish	INMETRO-certifikat för elektrisk utrustning för medicinskt bruk				
Turkish	INMETRO Sertifikası Tıbbi Elektrikli Ekipman				


Russia and Eurasia

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	Conformity requirements for the Russian Federation	Indicates that the <i>medical device</i> complies with the requirements established in the Russian GOST-R standards	English	Conformity mark GOST-R
				Bulgarian	Маркировка за съответствие GOST-R
				Croatian	Oznaka sukladnosti GOST-R
				Chinese	合规标志 GOST-R
				Czech	Značka shody GOST-R
				Danish	Overensstemmelsesmærke GOST-R
				Dutch	Conformiteitsmerkteken GOST-R
				Estonian	Vastavusmärgis GOST-R
				Finnish	GOST-R-sertifiointi
				French	Marque de conformité GOST-R
				German	Konformitätszeichen GOST-R
				Greek	Σήμα συμμόρφωσης GOST-R
				Hungarian	GOST-R-megfelelőségi jelölés
				Italian	Marchio di conformità GOST-R
				Japanese	適合マーク GOST-R
				Latvian	GOST-R atbilstības marķējums
				Lithuanian	Atitikties ženklas GOST-R
				Norwegian	Samsvarsmerke GOST-R
				Polish	Znak zgodności GOST-R
				Portuguese	Marca de conformidade GOST-R
Romanian	Marcaj de conformitate GOST-R				
Slovak	Značka zhody GOST-R				
Slovenian	Oznaka skladnosti GOST-R				
Spanish	Marcado de conformidad con GOST-R				
Swedish	Märkning om överensstämmelse, GOST-R				
Turkish	GOST-R uygunluk işareti				


Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	Technical Regulations of the Eurasian Economic Union (EAEU commission No. 711)	Indicates that the medical device meets the regulations and standards for customs clearance of the Eurasian Customs Union	English	Conformity mark to technical regulations of the Eurasian Customs Union
				Bulgarian	Маркировка за съответствие с техническите разпоредби на Евразийския митнически съюз
				Croatian	Oznaka sukladnosti s tehničkim propisima Euroazijske carinske unije
				Chinese	欧亚关税同盟技术法规的合规标志
				Czech	Značka shody s technickými předpisy Euroasijské celní unie
				Danish	Overensstemmelsesmærke i henhold til de tekniske forskrifter i Den Eurasiske Toldunion
				Dutch	Markering van conformiteit met technische voorschriften van de Euraziatische douane-unie
				Estonian	Vastavusmärgis, mis tähistab vastavust Euraasia Tolliliidu tehnilistele eeskirjadele
				Finnish	Euraasian tulliliiton teknisten määräysten yhdenmukaisuusmerkki
				French	Marque de conformité aux réglementations techniques de l'Union douanière eurasiatique
				German	Konformitätszeichen mit technischen Vorschriften der Eurasischen Zollunion.
				Greek	Σήμα συμμόρφωσης προς τους τεχνικούς κανονισμούς της Ευρασιατικής Τελωνειακής Ένωσης
				Hungarian	Az Eurázsiai Vámunió műszaki előírásainak való megfeleléség jelölése


				Italian	Marchio di conformità ai regolamenti tecnici dell'Unione doganale eurasiatica
				Japanese	ユーラシア関税同盟の技術規則への適合マーク
				Latvian	Zīme, kas apliecina atbilstību Eirāzijas muitas savienības tehniskajiem noteikumiem
				Lithuanian	Atitikties Eurazijos muitų sąjungos techniniams reglamentams ženklas
				Norwegian	Samsvarsmerke i henhold til tekniske forskrifter i den eurasiske tollunionen
				Polish	Znak zgodności z przepisami technicznymi Euroazjatyckiej Unii Celnej
				Portuguese	Marca de conformidade com os regulamentos técnicos da União Aduaneira Eurasiática
				Romanian	Marcaj de conformitate cu reglementările tehnice ale Uniunii Vamale Eurasiatice
				Slovak	Značka zhody s technickými predpismi Eurázijskej colnej únie
				Slovenian	Oznaka skladnosti s tehničnimi predpisi Evrazijske carinske unije
				Spanish	Marca de conformidad con los reglamentos técnicos de la Unión Aduanera de Eurasia
				Swedish	Överensstämmelsemärkning enligt eurasiska tullunionens tekniska föreskrifter
				Turkish	Avrasya Gümrük Birliği'nin teknik talimatlarına uygunluk işareti.

Gulf Countries – Bahrain, Kuwait, Oman, Qatar, Saudi Arabia and United Arab Emirates



Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 G. H. C.	N/A	N/A	Indicates that the product is in compliance with medical device requirements in the Gulf Cooperation Council countries	English	Gulf Health Council
				Bulgarian	Съвет по здравеопазване в Персийския залив
				Croatian	Zdravstveno vijeće GCC-a
				Chinese	海湾国家健康委员会
				Czech	Rada pro zdraví v Perském zálivu
				Danish	Golfstateres sundhedsråd
				Dutch	Gulf Health Council
				Estonian	Pärsia lahe tervisenõukogu
				Finnish	GHC, GCC- maiden terveysterveysto
				French	Conseil de la santé du Golfe (GHC)
				German	Golfstaaten-Gesundheitsrat
				Greek	Συμβούλιο Υγείας του Κόλπου
				Hungarian	Öböl Menti Egészségügyi Tanács
				Italian	Gulf Health Council
				Japanese	湾岸保健協議会 (Gulf Health Council)
				Latvian	Persijas līča Veselības padome
				Lithuanian	Persijos įlankos sveikatos taryba
				Norwegian	Gulf Health Council
				Polish	Gulf Health Council
				Portuguese	Conselho de Saúde do Golfo
				Romanian	Consiliul de Sănătate din Golf
Slovak	Rada pre zdravie arabských štátov v Perzskom zálive				
Slovenian	Svet za zdravje v zalivskih državah				
Spanish	Consejo de Cooperación del Golfo				
Swedish	Gulf Health Council				
Turkish	Körfez Sağlık Konseyi				

Miscellaneous

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Indicates a number used to identify trade items at various packaging levels	English	Global Trade Item Number
				Bulgarian	Глобален номер на търговската единица (GTIN)
				Croatian	Globalni broj trgovačke jedinice
				Chinese	全球贸易项目代码
				Czech	Číslo globální obchodní položky (GTIN)
				Danish	Globalt varenummer (GTIN)
				Dutch	Wereldwijd handelsartikelnummer (GTIN)
				Estonian	Globaalne kaubaartikli number
				Finnish	GTIN-koodi
				French	Code article international (GTIN)
				German	Globale Artikelnummer (GTIN)
				Greek	Διεθνής Κωδικός Μονάδων Εμπορίας (Global Trade Item Number)
				Hungarian	Globális kereskedelmi áruazonosító szám (GTIN)
				Italian	Codice articolo commerciale globale (GTIN)
				Japanese	グローバル・トレード・アイテム・ナンバー
				Latvian	Globālais tirdzniecības preces numurs
				Lithuanian	Pasaulinis prekės identifikavimo numeris
				Norwegian	Globalt handelsvarenummer (GTIN)
				Polish	Globalny numer identyfikacyjny jednostki handlowej (GTIN)
				Portuguese	Número de identificação de comércio global
Romanian	Număr global de articol comercial (GTIN)				
Slovak	Číslo globálnej obchodnej položky (GTIN)				
Slovenian	Globalna trgovinska številka izdelka				
Spanish	Número de artículo comercial internacional				
Swedish	Globalt artikelnummer (GTIN)				
Turkish	Küresel Ticaret Ürün Numarası				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
	N/A	N/A	Standard diameter symbol	English	Tube diameter
				Bulgarian	Диаметър на тръбата
				Croatian	Promjer cijevi
				Chinese	管径
				Czech	Průměr trubice
				Danish	Slangediameter
				Dutch	Diameter van de slang
				Estonian	Toru/vooliku diameeter
				Finnish	Putken halkaisija
				French	Diamètre du tube
				German	Tubusdurchmesser
				Greek	Διάμετρος σωλήνα
				Hungarian	Csőátmérő
				Italian	Diametro del tubo
				Japanese	チューブ半径
				Latvian	Stobriņa diametrs
				Lithuanian	Vamzdelio skersmuo
				Norwegian	Tubediameter
				Polish	Średnica rury
				Portuguese	Diâmetro do tubo
Romanian	Diametrul tubului				
Slovak	Priemer hadičky				
Slovenian	Premer cevke				
Spanish	Diámetro del tubo				
Swedish	Slangdiameter				
Turkish	Boru çapı				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
$\frac{\text{Ag}}{\text{AgCl}}$	N/A	N/A	Indicates that the product contains silver and silver chloride	English	Silver/silver chloride
				Bulgarian	Сребро/сребърен хлорид
				Croatian	Srebro/srebrov klorid
				Chinese	银/氯化银
				Czech	Stříbro/chlorid stříbrný
				Danish	Sølv/sølvklorid
				Dutch	Zilver/zilverchloride
				Estonian	Hõbe/hõbekloriid
				Finnish	Hopea/hopeakloridi
				French	Argent/Chlorure d'argent
				German	Silber/Silberchlorid
				Greek	Άργυρος/χλωριούχος άργυρος
				Hungarian	Ezüst/ezüst-klorid
				Italian	Argento/cloruro d'argento
				Japanese	銀/塩化銀
				Latvian	Sudrabs/sudraba hlorīds
				Lithuanian	Sidabras / sidabro chloridas
				Norwegian	Sølv/sølvklorid
				Polish	Srebro/chlorek srebra
				Portuguese	Prata/cloreto de prata
Romanian	Argint/clorură de argint				
Slovak	Striebro/chlorid strieborný				
Slovenian	Srebro/srebrov klorid				
Spanish	Plata / cloruro de plata				
Swedish	Silver/silverklorid				
Turkish	Gümüş/gümüş klorür				

Symbol	Symbol reference number	Standard – Title and Designation	Description	Language	Title
 	N/A	N/A	Indicates that the <i>medical device</i> does not contain pvc	English	Product is PVC free
				Bulgarian	Продуктът не съдържа PVC
				Croatian	Proizvod ne sadrži PVC
				Chinese	产品不含 PVC
				Czech	Výrobek neobsahuje PVC
				Danish	Produktet er PVC-frit
				Dutch	Product is pvc-vrij
				Estonian	Toode ei sisalda PVC-d
				Finnish	Tuote ei sisällä PVC:tä
				French	Le produit ne contient pas de PVC
				German	Produkt ist frei von PVC
				Greek	Το προϊόν δεν περιέχει PVC
				Hungarian	A termék PVC-mentes
				Italian	Il prodotto non contiene PVC
				Japanese	製品には PVC は使用されていません
				Latvian	Izstrādājums nesatur PVC
				Lithuanian	Gaminyje nėra PVC
				Norwegian	Produktet er PVC-fritt
				Polish	Produkt nie zawiera PCW
				Portuguese	O produto não contém PVC
Romanian	Produsul nu conține PVC				
Slovak	Výrobok neobsahuje PVC				
Slovenian	Izdelek ne vsebuje PVC				
Spanish	El producto no contiene PVC				
Swedish	Produkten är PVC-fri				
Turkish	Ürün PVC içermez				